

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food - Oils and Fats for Babyfood

Detection technique: LC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 10

QTI list nr.: 11

Version nr.: 3

Date issued: 15-08-2022

Valid till: N/A, current version



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Legend for excluded commodities

-

Component/Summation and CASnr	Reporting Limit mg/kg	Excluded
Q 2,3,5-Trimethacarb / 3,4,5-Trimethacarb (2655-15-4 / 2686-99-9)	0,001	
2,4,5-T (93-76-5)	0,010	
2,4,5-TP, Fenoprop, Fenoprop-butyl (93-72-1)	0,005	
2,4-D (94-75-7)	0,005	
Q 3-Chloroaniline (108-42-9)	0,005	
Q 3-Hydroxycarbofuran (16655-82-6)	0,001	
Q 4-Chloro-2-methylaniline (95-69-2)	0,001	
Q 4-CPA (122-88-3)	0,010	
4-Nitrophenol (100-02-7)	0,005	
Q Acephate (30560-19-1)	0,001	
Q Acetamiprid (135410-20-7)	0,001	
Q Acetamiprid, N-desmethyl- (190604-92-3)	0,001	
Acetochlor (34256-82-1)	0,010	
Acibenzolar Acid, CGA 210007 (35272-27-6)	0,005	
Q Acibenzolar-S-Methyl (135158-54-2)	0,001	
Q Aclonifen (74070-46-5)	0,005	
Q Albendazole (54965-21-8)	0,001	
Q Aldicarb (116-06-3)	0,001	
Q Aldicarb-Sulfone, Aldoxycarb (1646-88-4)	0,001	
Aldicarb-Sulfoxide (1646-87-3)	0,010	
Ametoctradin (865318-97-4)	0,001	
Q Amidosulfuron (120923-37-7)	0,001	
Q Aminocarb (2032-59-9)	0,001	
Q Aminopyralid (150114-71-9)	0,010	
Q Amisulbrom (348635-87-0)	0,002	
Q Anilofos (64249-01-0)	0,001	
Q Aspon (3244-90-4)	0,001	
Q Asulam (3337-71-1)	0,001	
Q Atraton (1610-17-9)	0,001	
Q Atrazine (1912-24-9)	0,001	
Q Atrazine-desethyl (6190-65-4)	0,001	
Q Azimsulfuron (120162-55-2)	0,001	
Q Azinphos-ethyl (2642-71-9)	0,001	
Q Azinphos-methyl (86-50-0)	0,001	
Q Aziprotryne (4658-28-0)	0,001	
Azoxystrobin (131860-33-8)	0,005	
BAC-10 (965-32-2)	0,001	
BAC-12 (139-07-1)	0,005	
BAC-14 (139-08-2)	0,005	
BAC-16 (122-18-9)	0,001	
Barban (101-27-9)	0,010	
Q Benazolin (3813-05-6)	0,001	
Q Benazolin-Ethyl (25059-80-7)	0,001	
Q Bendiocarb (22781-23-3)	0,001	
Q Bensulfuron-Methyl (83055-99-6)	0,001	
Bentazone (25057-89-0)	0,005	
Q Benthiavalicarb-Isopropyl (177406-68-7)	0,001	
Q Benthiazole (21564-17-0)	0,001	
Q Benzobicyclon (156963-66-5)	0,001	

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food - Oils and Fats for Babyfood

Detection technique: LC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 10

QTI list nr.: 11

Version nr.: 3

Date issued: 15-08-2022

Valid till: N/A, current version



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Legend for excluded
commodities

-

Component/Summation and CASnr	Reporting Limit mg/kg	Excluded
Q Benzovindiflupyr (Solatenol) (1072957-71-1)	0,001	
Bifenazate (149877-41-8)	0,001	
Bifenazate-diazene (149878-40-0)	0,001	
Q Bispyribac (125401-75-4)	0,001	
Q Bixafen (581809-46-3)	0,001	
Q Boscalid (188425-85-6)	0,010	
Brodifacoum (56073-10-0)	0,005	
Q Bromacil (314-40-9)	0,001	
Q Bromoxynil (1689-84-5)	0,001	
Bromuconazole a (116255-48-2)	0,001	
Bromuconazole b (116255-48-2)	0,001	
Q Bupirimate (41483-43-6)	0,001	
Q Buprofezin (69327-76-0)	0,001	
Q Butachlor (23184-66-9)	0,001	
Butafenacil (134605-64-4)	0,010	
Q Butocarboxim (34681-10-2)	0,001	
Q Butocarboxim-Sulfoxide (34681-24-8)	0,005	
Q Butoxycarboxim (34681-23-7)	0,001	
Q Buturon (3766-60-7)	0,001	
Q Butylate (2008-41-5)	0,001	
Q Cadusafos (95465-99-9)	0,001	
Q Carbaryl (63-25-2)	0,005	
Q Carbazim (10605-21-7)	0,001	
Q Carbetamide (16118-49-3)	0,001	
Carbofuran (1563-66-2)	0,001	
Q Carbofuran-3-keto (16709-30-1)	0,001	
Carbophenothion (786-19-6)	0,005	
Q Carboxin (5234-68-4)	0,010	
Q Carfentrazone-Ethyl (128639-02-1)	0,001	
Q Carpropamid (104030-54-8)	0,001	
Chinomethionat (2439-01-2)	0,001	
Q Chlorantraniliprole (Rynaxapyr) (500008-45-7)	0,001	
Chlorbufam (1967-16-4)	0,005	
Q Chlordimeform (6164-98-3)	0,001	
Q Chlorfenvinphos (470-90-6)	0,001	
Q Chlorfluazuron (71422-67-8)	0,001	
Chlorflurenol-Methyl ester (2536-31-4)	0,010	
Q Chloridazon (1698-60-8)	0,001	
Q Chlorimuron-Ethyl (90982-32-4)	0,001	
Q Chlormequat (7003-89-6)	0,005	
Q Chlorobenzuron (57160-47-1)	0,001	
Q Chloroxuron (1982-47-4)	0,001	
Q Chlorpyrifos (2921-88-2)	0,001	
Q Chlorpyrifos-Methyl (5598-13-0)	0,001	
Q Chlorsulfuron (64902-72-3)	0,001	
Q Chlortoluron (15545-48-9)	0,001	
Q Chromafenozide (143807-66-3)	0,001	
Q Cinidon-Ethyl (142891-20-1)	0,010	
Q Cinosulfuron (94593-91-6)	0,001	

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food - Oils and Fats for Babyfood

Detection technique: LC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 10

QTI list nr.: 11

Version nr.: 3

Date issued: 15-08-2022

Valid till: N/A, current version



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Legend for excluded commodities

-

Component/Summation and CASnr	Reporting Limit mg/kg	Excluded
Clethodim (99129-21-2)	0,001	
Q Clethodim-Sulfone (111031-17-5)	0,001	
Clodinafop (114420-56-3)	0,020	
Q Clodinafop-Propargyl (105512-06-9)	0,001	
Clofentezine (74115-24-5)	0,001	
Q Clomazone (81777-89-1)	0,001	
Q Clopyralid (1702-17-6)	0,005	
Q Clothianidin (210880-92-5)	0,001	
Q Coumaphos (56-72-4)	0,001	
Coumatetralyl (5836-29-3)	0,001	
Q Crimidine (535-89-7)	0,001	
Q Crotoxyphos (7700-17-6)	0,001	
Crufomate (299-86-5)	0,001	
Cyantraniliprole (736994-63-1)	0,001	
Q Cyazofamid (120116-88-3)	0,001	
Q Cyclanilide (113136-77-9)	0,010	
Cycloheximide (66-81-9)	0,005	
Q Cyclosulfamuron (136849-15-5)	0,001	
Cycloxydim (101205-02-1)	0,001	
Q Cycluron, OMU (2163-69-1)	0,005	
Q Cyenopyrafen (560121-52-0)	0,001	
Q Cyflufenamid: sum of cyflufenamid (Z-isomer) and its E-isomer (180409-60-3)	0,001	
Q Cyflumetofen (400882-07-7)	0,001	
Q Cymoxanil (57966-95-7)	0,001	
Q Cyproconazole A (94361-06-5)	0,001	
Q Cyproconazole B (94361-06-5)	0,001	
Q Cyprodinil (121552-61-2)	0,005	
Cyromazine (66215-27-8)	0,005	
Q Daimuron (42609-52-9)	0,001	
Q Demeton-O (298-03-3)	0,001	
Q Demeton-S (126-75-0)	0,001	
Demeton-s-methyl (919-86-8)	0,001	
Q Demeton-S-Methyl-sulfon (17040-19-6)	0,001	
Desmedipham (13684-56-5)	0,001	
Di-Allate (2303-16-4)	0,050	
Q Diazinon (333-41-5)	0,001	
Q Dichlorimid (37764-25-3)	0,010	
Dichlorprop, 2,4-DP (120-36-5)	0,001	
Dichlorvos (62-73-7)	0,001	
Q Diclobutrazol (75736-33-3)	0,001	
Diclofop (40843-25-2)	0,010	
Q Diclofop-Methyl (51338-27-3)	0,010	
Q Dicrotophos (141-66-2)	0,001	
Diethatyl-Ethyl (38727-55-8)	0,001	
Q Difenoconazole (119446-68-3)	0,001	
Q Difenoxuron (14214-32-5)	0,001	
Q Difenzoquat methyl sulfate (43222-48-6)	0,005	
Q Diflubenzuron (35367-38-5)	0,001	
Q Diflufenican (83164-33-4)	0,001	

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food - Oils and Fats for Babyfood

Detection technique: LC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 10

QTI list nr.: 11

Version nr.: 3

Date issued: 15-08-2022

Valid till: N/A, current version



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Legend for excluded
commodities

-

Component/Summation and CASnr	Reporting Limit mg/kg	Excluded
Q Dimefox (115-26-4)	0,001	
Q Dimefuron (34205-21-5)	0,001	
Dimethachlor (50563-36-5)	0,001	
Q Dimethenamid (87674-68-8)	0,001	
Q Dimethoate (60-51-5)	0,001	
Q Dimethomorph (110488-70-5)	0,001	
Dimethylphenylformamide, 2,4-, DMF (60397-77-5)	0,010	
Dimethylphenyl-N-Methylformamidine, N-2,4-, DMPF (33089-74-6)	0,050	
Q Dimetilan (644-64-4)	0,001	
Q Dimoxystrobin (149961-52-4)	0,001	
Q Diniconazole (83657-24-3)	0,001	
Q Dinotefuran (165252-70-0)	0,005	
Q Dinoterb, DNTBP (1420-07-1)	0,010	
Q Dioxabenzofos (Salithion) (3811-49-2)	0,005	
Q Dioxacarb (6988-21-2)	0,001	
Diphenylamine (122-39-4)	0,020	
Q Disulfoton-sulfone (2497-06-5)	0,001	
Q Disulfoton-sulfoxide (2497-07-6)	0,001	
Q Ditalimfos (5131-24-8)	0,001	
Q Dithiopyr (97886-45-8)	0,001	
Q Diuron (330-54-1)	0,001	
Q DMS (3984-14-3)	0,005	
Q DNOC (534-52-1)	0,005	
Q Dodemorph (1593-77-7)	0,001	
Q EPN (2104-64-5)	0,010	
Q Epoxiconazole (133855-98-8)	0,001	
Q EPTC (759-94-4)	0,005	
Q Esprocarb (85785-20-2)	0,001	
Q Etaconazole (60207-93-4)	0,010	
Q Ethaboxam (162650-77-3)	0,001	
Q Ethametsulfuron-Methyl (97780-06-8)	0,001	
Q Ethiofencarb-Sulfone (53380-23-7)	0,005	
Q Ethiofencarb-Sulfoxide (53380-22-6)	0,001	
Q Ethiprole (181587-01-9)	0,005	
Q Ethirimol (23947-60-6)	0,005	
Q Ethoprophos (13194-48-4)	0,001	
Q Ethoxysulfuron (126801-58-9)	0,001	
Q Etrimfos (38260-54-7)	0,001	
Q Famoxadone (131807-57-3)	0,001	
Q Famphur, Famophos (52-85-7)	0,001	
Q Fenamidone (161326-34-7)	0,001	
Fenamiphos (22224-92-6)	0,001	
Q Fenamiphos-Sulfoxide (31972-43-7)	0,001	
Q Fenarimol (60168-88-9)	0,001	
Q Fenazaquin (120928-09-8)	0,001	
Q Fenbuconazole (114369-43-6)	0,001	
Q Fenchlorazole-Ethyl (103112-35-2)	0,001	
Q Fenchlorphos oxon (3983-45-7)	0,001	
Q Fenhexamid (126833-17-8)	0,005	

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food - Oils and Fats for Babyfood

Detection technique: LC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 10

QTI list nr.: 11

Version nr.: 3

Date issued: 15-08-2022

Valid till: N/A, current version



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Legend for excluded commodities

-

Component/Summation and CASnr	Reporting Limit mg/kg	Excluded
Q Fenobucarb (3766-81-2)	0,001	
Q Fenoxanil (115852-48-7)	0,001	
Q Fenoxaprop and Fenoxaprop-P (95617-09-7 / 113158-40-0)	0,005	
Q Fenoxaprop-Ethyl and Fenoxaprop-P-Ethyl (66441-23-4 / 71283-80-2)	0,001	
Q Fenoxycarb (72490-01-8)	0,001	
Q Fenpiclonil (74738-17-3)	0,001	
Q Fenpropidin (67306-00-7)	0,001	
Q Fenpropimorph (67564-91-4)	0,001	
Q Fenpyrazamine (473798-59-3)	0,001	
Q Fenpyroximate (134098-61-6/111812-58-9)	0,001	
Q Fensulfothion (115-90-2)	0,001	
Q Fensulfothion-Oxon (6552-21-2)	0,001	
Q Fensulfothion-oxon-sulfon (6132-17-8)	0,001	
Q Fensulfothion-sulfone (14255-72-2)	0,001	
Fenthion (55-38-9)	0,005	
Q Fenthion-Oxonsulfone (14086-35-2)	0,001	
Q Fenthion-Oxonsulfoxide (6552-13-2)	0,001	
Q Fenthion-Sulfone (3761-42-0)	0,001	
Q Fenuron (101-42-8)	0,001	
Q Ferimzone (sum of (E),(Z) types) (89269-64-7)	0,005	
Fipronil (120068-37-3)	0,005	
Fipronil-Sulfone (120068-36-2)	0,050	
Q Flamprop-M-isopropyl (57973-67-8)	0,005	
Q Flonicamid (158062-67-0)	0,001	
Flonicamid-Metabolite TFNG (207502-65-6)	0,005	
Q Florasulam (145701-23-1)	0,001	
Q Fluacrypyrim (229977-93-9)	0,005	
Q Fluazifop and Fluazifop-P (69335-91-7 / 83066-88-0)	0,001	
Q Fluazifop-butyl and Fluazifop-P-butyl (69806-50-4 / 79241-46-6)	0,010	
Fluazinam (79622-59-6)	0,001	
Fluazuron (86811-58-7)	0,001	
Flubendiamide (272451-65-7)	0,010	
Q Flucetosulfuron (412928-75-7)	0,001	
Q Flucycloxuron (113036-88-7)	0,001	
Fludioxonil (131341-86-1)	0,001	
Q Flufenacet (142459-58-3)	0,001	
Q Flufenoxuron (101463-69-8)	0,001	
Q Flumetsulam (98967-40-9)	0,001	
Q Flumiclorac-pentyl (87546-18-7)	0,001	
Flumioxazin (103361-09-7)	0,005	
Q Fluometuron (2164-17-2)	0,001	
Q Fluopicolide (239110-15-7)	0,010	
Q Fluopyram (658066-35-4)	0,001	
Q Fluoroglycofene-Ethyl (77501-90-7)	0,005	
Q Fluoxastrobin (361377-29-9)	0,001	
Q Flupyradifurone (951659-40-8)	0,001	
Q Flupyrasulfuron-methyl (144740-53-4)	0,001	
Q Fluquinconazole (136426-54-5)	0,001	
Q Fluralaner (864731-61-3)	0,001	

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food - Oils and Fats for Babyfood

Detection technique: LC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 10

QTI list nr.: 11

Version nr.: 3

Date issued: 15-08-2022

Valid till: N/A, current version



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Legend for excluded
commodities

-

Component/Summation and CASnr	Reporting Limit mg/kg	Excluded
Q Fluridone (59756-60-4)	0,001	
Q Flurochloridone (61213-25-0)	0,001	
Q Fluroxypyr (69377-81-7)	0,005	
Fluroxypyr-Meptyl (81406-37-3)	0,001	
Q Flurtamone (96525-23-4)	0,005	
Q Flusilazole (85509-19-9)	0,001	
Q Flusulfamide (106917-52-6)	0,001	
Q Fluthiacet-Methyl (117337-19-6)	0,010	
Q Flutriafol (76674-21-0)	0,001	
Q Fluxapyroxad (907204-31-3)	0,001	
Fomesafen (72178-02-0)	0,005	
Q Fonofos (944-22-9)	0,001	
Q Foramsulfuron (173159-57-4)	0,001	
Q Forchlorfenuron (68157-60-8)	0,001	
Q Formetanate (22259-30-9)	0,001	
Q Formothion (2540-82-1)	0,005	
Q Fosthiazate (98886-44-3)	0,001	
Q Fuberidazole (3878-19-1)	0,001	
Q Furalaxyl (57646-30-7)	0,005	
Furathiocarb (65907-30-4)	0,001	
Furmecyclox (60568-05-0)	0,001	
Q Halofenozide (112226-61-6)	0,001	
Q Halosulfuron-methyl (100784-20-1)	0,001	
Q Haloxyfop (69806-34-4)	0,001	
Q Haloxyfop-2-ethoxyethyl (87237-48-7)	0,001	
Haloxyfop-methyl and Haloxyfop-p-Methyl (69806-40-2 / 72619-32-0)	0,001	
Q Heptenophos (23560-59-0)	0,005	
Q Hexaconazole (79983-71-4)	0,001	
Hexaflumuron (86479-06-3)	0,005	
Q Hexazinone (51235-04-2)	0,001	
Q Hexythiazox (78587-05-0)	0,001	
Q Hydramethylnon (67485-29-4)	0,020	
Hymexazol, Hydroxyisoxazol (10004-44-1)	0,005	
Q Imazalil (35554-44-0)	0,010	
Q Imazamethabenz-Methyl (81405-85-8)	0,001	
Q Imazamox (114311-32-9)	0,010	
Q Imazapic (104098-48-8)	0,005	
Q Imazapyr (81334-34-1)	0,005	
Q Imazaquin (81335-37-7)	0,001	
Q Imazethapyr (81335-77-5)	0,001	
Q Imazosulfuron (122548-33-8)	0,001	
Q Imibenconazole (86598-92-7)	0,001	
Q Imidacloprid (138261-41-3)	0,001	
Indole-3-Acetic Acid, IAA (87-51-4)	0,005	
Q Indole-3-Butyric Acid, IBA (133-32-4)	0,001	
Q Indoxacarb (144171-61-9)	0,002	
Q Ioxynil (1689-83-4)	0,001	
Q Ipconazole (125225-28-7)	0,001	
Q Iprobenfos (26087-47-8)	0,001	

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food - Oils and Fats for Babyfood

Detection technique: LC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 10

QTI list nr.: 11

Version nr.: 3

Date issued: 15-08-2022

Valid till: N/A, current version



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Legend for excluded
commodities

-

Component/Summation and CASnr	Reporting Limit mg/kg	Excluded
Q Iprodione (36734-19-7)	0,001	
Q Iprovalicarb (140923-17-7)	0,001	
Q Isazofos (42509-80-8)	0,001	
Q Isofenphos (25311-71-1)	0,001	
Q Isofenphos-Methyl (99675-03-3)	0,005	
Q Isofetamid (875915-78-9)	0,001	
Q Isomethiozin (57052-04-7)	0,001	
Q Isoprocab (2631-40-5)	0,001	
Isopropalin (33820-53-0)	0,001	
Q Isoprothiolane (50512-35-1)	0,001	
Q Isoproturon (34123-59-6)	0,001	
Q Isopyrazam (881685-58-1)	0,001	
Q Isoxaben (82558-50-7)	0,001	
Q Isoxaflutole (141112-29-0)	0,001	
Karanjin (521-88-0)	0,005	
Q Karbutilate (4849-32-5)	0,001	
Q Kresoxim-methyl (143390-89-0)	0,001	
Lenacil (2164-08-1)	0,005	
Q Linuron (330-55-2)	0,001	
Q Lufenuron (103055-07-8)	0,001	
Q Malaoxon (1634-78-2)	0,001	
Malathion (121-75-5)	0,005	
Q Mandestrobin (173662-97-0)	0,001	
Q Mandipropamid (374726-62-2)	0,001	
MCPA (94-74-6)	0,005	
MCP, Mechlorprop, Mecoprop (7085-19-0)	0,001	
Q Mecarbam (2595-54-2)	0,001	
Q Mefenpyr-diethyl (135590-91-9)	0,001	
Q Mepanipyrim (110235-47-7)	0,005	
Q Mepiquat (15302-91-7)	0,001	
Q Mesosulfuron-Methyl (208465-21-8)	0,001	
Mesotrione (104206-82-8)	0,001	
Q Metaflumizone (139968-49-3)	0,001	
Q Metalaxyl and Metalaxyl-M (57837-19-1 / 70630-17-0)	0,001	
Metalaxyl free acid (87764-37-2)	0,005	
Q Metamitron (41394-05-2)	0,001	
Q Metazachlor (67129-08-2)	0,001	
Q Metconazole (125116-23-6)	0,001	
Q Methabenzthiazuron (18691-97-9)	0,001	
Q Methamidophos (10265-92-6)	0,005	
Methidathion (950-37-8)	0,001	
Q Methiocarb (2032-65-7)	0,001	
Q Methiocarb-sulfon (2179-25-1)	0,001	
Q Methiocarb-sulfoxide (2635-10-1)	0,001	
Q Metholachlor (51218-45-2)	0,001	
Methomyl (16752-77-5)	0,010	
Q Methoprotetryne (841-06-5)	0,001	
Q Methoxyfenozide (161050-58-4)	0,001	
Q Metobromuron (3060-89-7)	0,001	

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food - Oils and Fats for Babyfood

Detection technique: LC-MS/MS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 10

QTI list nr.: 11

Version nr.: 3

Date issued: 15-08-2022

Valid till: N/A, current version



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Legend for excluded
commodities

-

Component/Summation and CASnr	Reporting Limit mg/kg	Excluded
Metominostrobin E (133408-50-1)	0,005	
Q Metosulam (139528-85-1)	0,001	
Q Metoxuron (19937-59-8)	0,001	
Q Metrafenone (220899-03-6)	0,001	
Metribuzin (21087-64-9)	0,001	
Q Metsulfuron-Methyl (74223-64-6)	0,001	
Q Mexacarbate (315-18-4)	0,001	
MGK-264 (113-48-4)	0,001	
Q Molinate (2212-67-1)	0,001	
Monocrotophos (6923-22-4)	0,001	
Q Monolinuron (1746-81-2)	0,001	
Q Monuron (150-68-5)	0,001	
Q Napropamide (15299-99-7)	0,001	
Q Neburon (555-37-3)	0,001	
Q Nicosulfuron (111991-09-4)	0,001	
Q Nitenpyram (150824-47-8)	0,001	
Q Nitralin (4726-14-1)	0,010	
Q Novaluron (116714-46-6)	0,001	
Q Ofurace (58810-48-3)	0,001	
Q Omethoate (1113-02-6)	0,001	
Q Orbencarb (34622-58-7)	0,001	
Q Orthosulfamuron (213464-77-8)	0,020	
Q Oryzalin (19044-88-3)	0,001	
Q Oxadiargyl (39807-15-3)	0,005	
Q Oxadiazon (19666-30-9)	0,001	
Q Oxadixyl (77732-09-3)	0,001	
Oxamyl (23135-22-0)	0,001	
Q Oxamyl-Oxime (30558-43-1)	0,001	
Q Oxasulfuron (144651-06-9)	0,001	
Q Oxathiapiprolin (1003318-67-9)	0,001	
Q Oxaziclonofone (153197-14-9)	0,001	
Q Oxycarboxin (5259-88-1)	0,001	
Q Oxydemeton-Methyl (301-12-2)	0,001	
Q Paclobutrazol (76738-62-0)	0,001	
Q Pebulate (1114-71-2)	0,001	
Q Penconazole (66246-88-6)	0,001	
Q Pencycuron (66063-05-6)	0,001	
Q Penoxsulam (219714-96-2)	0,001	
Pentachlorophenol (87-86-5)	0,005	
Q Penthiopyrad (183675-82-3)	0,001	
Q Pethoxamid (106700-29-2)	0,001	
Phenmedipham (13684-63-4)	0,001	
Q Phenthoate (2597-03-7)	0,010	
Phorate-Oxon (2600-69-3)	0,002	
Q Phorate-Oxonsulfone (2588-06-9)	0,001	
Q Phorate-Oxonsulfoxide (2588-05-8)	0,001	
Q Phorate-sulfone (2588-04-7)	0,005	
Q Phorate-sulfoxide (2588-03-6)	0,001	
Q Phosmet (732-11-6)	0,001	

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food - Oils and Fats for Babyfood

Detection technique: LC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 10

QTI list nr.: 11

Version nr.: 3

Date issued: 15-08-2022

Valid till: N/A, current version



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Legend for excluded
commodities

-

Component/Summation and CASnr	Reporting Limit mg/kg	Excluded
Q Phosphamidon (13171-21-6)	0,001	
Q Phoxim (14816-18-3)	0,001	
Picloram (1918-02-1)	0,020	
Q Picolinafen (137641-05-5)	0,005	
Q Picoxystrobin (117428-22-5)	0,001	
Q Pinoxaden (243973-20-8)	0,001	
Q Piperonyl-butoxide (51-03-6)	0,001	
Q Pirimicarb (23103-98-2)	0,001	
Q Pirimicarb, Desmethyl- (30614-22-3)	0,001	
Pirimicarb, Desmethylformamido- (27218-04-8)	0,001	
Q Pirimiphos-Ethyl (23505-41-1)	0,001	
Q Pirimiphos-Methyl (29232-93-7)	0,001	
Q Pirimiphos-methyl-N-desethyl (67018-59-1)	0,001	
Pretilachlor (51218-49-6)	0,050	
Q Primisulfuron-Methyl (86209-51-0)	0,001	
Probenazole (27605-76-1)	0,005	
Q Prochloraz (67747-09-5)	0,001	
Q Prochloraz metabolite BTS44595 (139520-94-8)	0,001	
Q Prochloraz metabolite BTS44596 (139542-32-8)	0,001	
Q Prodiamine (29091-21-2)	0,005	
Q Profenofos (41198-08-7)	0,001	
Q Promecarb (2631-37-0)	0,001	
Q Prometryn (7287-19-6)	0,001	
Q Propachlor (1918-16-7)	0,001	
Propachlor oxalamic acid (70628-36-3)	0,005	
Q Propamocarb (24579-73-5)	0,001	
Q Propanil (709-98-8)	0,001	
Q Propaquizafop (111479-05-1)	0,010	
Q Propiconazole (60207-90-1)	0,001	
Q Propisochlor (86763-47-5)	0,001	
Propoxur (114-26-1)	0,001	
Propyzamide (23950-58-5)	0,001	
Proquinazid (189278-12-4)	0,001	
Q Prosulfuron (94125-34-5)	0,001	
Pymetrozine (123312-89-0)	0,005	
Q Pyracarbolid (24691-76-7)	0,001	
Q Pyraclonil (158353-15-2)	0,001	
Q Pyraclostrobin (175013-18-0)	0,001	
Q Pyrazolynate (58011-68-0)	0,001	
Q Pyrazophos (13457-18-6)	0,001	
Q Pyrazosulfuron-Ethyl (93697-74-6)	0,001	
Pyridalyl (179101-81-6)	0,005	
Q Pyridaphenthion (119-12-0)	0,001	
Q Pyridate (55512-33-9)	0,001	
Pyridate metabolite (40020-01-7)	0,001	
Q Pyrifenox (88283-41-4)	0,020	
Q Pyritalid (135186-78-6)	0,001	
Q Pyrimethanil (53112-28-0)	0,001	
Q Pyriminobac-methyl E (147411-69-6)	0,001	

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food - Oils and Fats for Babyfood

Detection technique: LC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 10

QTI list nr.: 11

Version nr.: 3

Date issued: 15-08-2022

Valid till: N/A, current version



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Legend for excluded
commodities

-

Component/Summation and CASnr	Reporting Limit mg/kg	Excluded
Q Pyriminobac-methyl Z (147411-70-9)	0,001	
Pyrimisulfan (221205-90-9)	0,002	
Q Pyriofenone (688046-61-9)	0,001	
Q Pyriproxyfen (95737-68-1)	0,001	
Q Pyroquilon, 4-Lilolidone (57369-32-1)	0,001	
Pyroxasulfone (447399-55-5)	0,005	
Q Pyroxsulam (422556-08-9)	0,001	
Q Quinalphos (13593-03-8)	0,001	
Quinclorac-methyl ester (84087-33-2)	0,010	
Quinmerac (90717-03-6)	0,001	
Q Quinoclamine, ACNQ (2797-51-5)	0,010	
Q Quinoxifen (124495-18-7)	0,001	
Q Quizalofop and Quizalofop-P (76578-12-6 / 94051-08-8)	0,002	
Q Quizalofop-Ethyl (76578-14-8)	0,001	
Resmethrin (10453-86-8 / 28434-01-7)	0,001	
Q Rimsulfuron (122931-48-0)	0,001	
Rotenone (83-79-4)	0,001	
Q Saflufenacil (372137-35-4)	0,001	
Q Sebuthylazine (7286-69-3)	0,001	
Q Sebuthylazine, Desethyl- (37019-18-4)	0,001	
Q Sedaxane (874967-67-6)	0,001	
Q Sethoxydim (74051-80-2)	0,001	
Q Siduron (1982-49-6)	0,001	
Q Simeconazole (149508-90-7)	0,001	
Q Simetryn (1014-70-6)	0,001	
Spinetoram A (187166-40-1)	0,001	
Spinetoram B (187166-40-1)	0,005	
Spinosyn A (168316-95-8)	0,001	
Spinosyn D (168316-95-8)	0,001	
Spirodiclofen (148477-71-8)	0,005	
Q Spiromesifen (283594-90-1)	0,001	
Spirotetramat (203313-25-1)	0,005	
Spiroxamine (sum) (118134-30-8)	0,001	
Sulcotrione (99105-77-8)	0,005	
Sulfentrazone (122836-35-5)	0,005	
Q Sulfometuron-Methyl (74222-97-2)	0,001	
Sulfosulfuron (141776-32-1)	0,001	
Q Sulfotep (3689-24-5)	0,001	
Sulfoxaflor (946578-00-3)	0,005	
Q Sulprofos (35400-43-2)	0,001	
Q Tebuconazole (107534-96-3)	0,001	
Q Tebufenozide (112410-23-8)	0,001	
Q Tebufenpyrad (119168-77-3)	0,001	
Tebutam (35256-85-0)	0,001	
Q Tebuthiuron (34014-18-1)	0,001	
Q Teflubenzuron (83121-18-0)	0,001	
Q TEPP (107-49-3)	0,001	
Q Tepraloxymid (149979-41-9)	0,001	
Q Terbufos-Sulfone (56070-16-7)	0,001	

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food - Oils and Fats for Babyfood

Detection technique: LC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 10

QTI list nr.: 11

Version nr.: 3

Date issued: 15-08-2022

Valid till: N/A, current version



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Legend for excluded
commodities

-

Component/Summation and CASnr	Reporting Limit mg/kg	Excluded
Terbufos-sulfoxide (10548-10-4)	0,001	
Q Terbumeton (33693-04-8)	0,001	
Q Terbutylazin-desethyl (30125-63-4)	0,001	
Q Terbutryn (886-50-0)	0,005	
Q Tetrachlorvinphos (22248-79-9)	0,001	
Q Tetraconazole (112281-77-3)	0,001	
Q Tetramethrin (7696-12-0)	0,001	
Q TFNA (158063-66-2)	0,005	
Q Thiabendazole (148-79-8)	0,001	
Q Thiabendazole, 5-Hydroxy- (948-71-0)	0,005	
Q Thiacloprid (111988-49-9)	0,001	
Q Thiamethoxam (153719-23-4)	0,005	
Q Thiazafuron (25366-23-8)	0,001	
Q Thiazopir (117718-60-2)	0,005	
Q Thidiazuron (51707-55-2)	0,001	
Q Thifensulfuron-methyl (79277-27-3)	0,005	
Q Thifluzamide (130000-40-7)	0,001	
Q Thiobencarb, Benthicarb (28249-77-6)	0,001	
Thiocyclam (31895-22-4)	0,005	
Q Thiofanox, Thiofanocarb (39196-18-4)	0,001	
Q Thiofanox-Sulfone (39184-59-3)	0,001	
Q Thiofanox-Sulfoxide (39184-27-5)	0,001	
Thiophanate-Ethyl (23564-06-9)	0,001	
Thiophanate-Methyl (23564-05-8)	0,001	
Tiadinil (223580-51-6)	0,010	
Q Tolclofos-Methyl (57018-04-9)	0,001	
Q Tolfenpyrad (129558-76-5)	0,001	
Topramezone (210631-68-8)	0,010	
Q Tralkoxydim (87820-88-0)	0,001	
Q Triadimefon (43121-43-3)	0,001	
Triadimenol (55219-65-3)	0,001	
Tri-allate (2303-17-5)	0,001	
Q Triasulfuron (82097-50-5)	0,001	
Tribenuron-Methyl (101200-48-0)	0,001	
Q Trichlorfon (52-68-6)	0,001	
Triclopyr (55335-06-3)	0,050	
Q Triclopyr-2-butoxyethyl ester (64700-56-7)	0,001	
Q Tricyclazole (41814-78-2)	0,001	
Q Tridemorph (24602-86-6)	0,005	
Q Trietazine (1912-26-1)	0,001	
Trifenylin (Fentin) (668-34-8)	0,002	
Q Trifloxystrobin (141517-21-7)	0,001	
Trifloxysulfuron (145099-21-4)	0,001	
Q Triflumizole (68694-11-1)	0,001	
Q Triflumuron (64628-44-0)	0,001	
Q Triflurosulfuron-Methyl (126535-15-7)	0,001	
Triforin (26644-46-2)	0,005	
Q Trinexapac-ethyl (95266-40-3)	0,001	
Q Triticonazole (131983-72-7)	0,001	

Current list of Pesticides within the scope of analytical protocol QTI-C-009:
Determination of the content of pesticides and additives
Equivalent to NEN-EN 15662
Scope: Food - Oils and Fats for Babyfood
Detection technique: LC-MSMS
Performance accredited scope in agreement with guideline SANTE/11321/2021
Annex to declaration of accreditation (scope of accreditation)
Normative document: EN ISO/IEC 17025:2017
Registration number RvA: L 678



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)
qti-services.com | qti@qti-services.com

Activity nr.: 10

QTI list nr.: 11

Version nr.: 3

Date issued: 15-08-2022

Valid till: N/A, current version

Legend for excluded
commodities

-

Component/Summation and CASnr	Reporting Limit mg/kg	Excluded
Q Tritosulfuron (142469-14-5)	0,001	
Q Valifenalate, Valiphenal (283159-90-0)	0,001	
Q Vamidothion (2275-23-2)	0,001	
Q Vernolate (1929-77-7)	0,010	
Q Zoxamide (156052-68-5)	0,001	

Mutations from previous version are highlighted in grey/dark