

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food and Feed - Fatty Acids

Detection technique: LC-MS/MS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 12

List nr.: 3

Version nr.: 9

Date issued: 28-05-2022

Valid till: 14-08-2022



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Component/Summation and CASnr	Reporting Limit mg/kg	Accredited (Y/N)
2,3,5-Trimethacarb / 3,4,5-Trimethacarb (2655-15-4 / 2686-99-9)	0,010	Y
2,4,5-T (93-76-5)	0,010	N
2,4,5-TP, Fenoprop, Fenoprop-butotyl (93-72-1)	0,005	N
2,4-D (94-75-7)	0,005	N
3-Chloroaniline (108-42-9)	0,005	Y
3-Hydroxycarbofuran (16655-82-6)	0,005	N
4-Nitrophenol (100-02-7)	0,005	N
Acephate (30560-19-1)	0,005	Y
Acetamiprid (135410-20-7)	0,010	N
Acetamiprid, N-desmethyl- (190604-92-3)	0,005	Y
Acibenzolar Acid, CGA 210007 (35272-27-6)	0,005	N
Acibenzolar-S-Methyl (135158-54-2)	0,010	Y
Albendazole (54965-21-8)	0,010	Y
Aldicarb (116-06-3)	0,010	Y
Aldicarb-Sulfone, Aldoxycarb (1646-88-4)	0,005	Y
Aldicarb-Sulfoxide (1646-87-3)	0,010	N
Amidosulfuron (120923-37-7)	0,010	Y
Aminocarb (2032-59-9)	0,005	Y
Anilofos (64249-01-0)	0,010	Y
Aspon (3244-90-4)	0,005	Y
Asulam (3337-71-1)	0,010	Y
Atraton (1610-17-9)	0,005	Y
Atrazine (1912-24-9)	0,005	Y
Atrazine-desethyl (6190-65-4)	0,005	Y
Azimsulfuron (120162-55-2)	0,020	Y
Azinphos-ethyl (2642-71-9)	0,010	Y
Azinphos-methyl (86-50-0)	0,010	Y
Azoxystrobin (131860-33-8)	0,005	N
BAC-12 (139-07-1)	0,005	N
BAC-14 (139-08-2)	0,005	N
BAC-16 (122-18-9)	0,005	N
Benazolin-Ethyl (25059-80-7)	0,005	Y
Bendiocarb (22781-23-3)	0,010	Y
Bensulfuron-Methyl (83055-99-6)	0,020	Y
Bentazone (25057-89-0)	0,005	N
Benthiavalicarb-Isopropyl (177406-68-7)	0,005	Y
Benthiazole (21564-17-0)	0,005	Y
Bifenazate (149877-41-8)	0,005	N
Boscalid (188425-85-6)	0,010	Y
Bromacil (314-40-9)	0,010	Y
Bromuconazole a (116255-48-2)	0,005	N
Bromuconazole b (116255-48-2)	0,005	N
Bupirimate (41483-43-6)	0,005	Y
Buprofezin (69327-76-0)	0,005	Y
Butafenacil (134605-64-4)	0,010	N
Butocarboxim (34681-10-2)	0,005	Y
Butocarboxim-Sulfoxide (34681-24-8)	0,005	Y
Butoxycarboxim (34681-23-7)	0,005	Y
Buturon (3766-60-7)	0,010	Y

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food and Feed - Fatty Acids

Detection technique: LC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 12

List nr.: 3

Version nr.: 9

Date issued: 28-05-2022

Valid till: 14-08-2022



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Component/Summation and CASnr	Reporting Limit mg/kg	Accredited (Y/N)
Butylate (2008-41-5)	0,005	Y
Carbaryl (63-25-2)	0,005	N
Carbendazim (10605-21-7)	0,005	N
Carbetamide (16118-49-3)	0,005	Y
Carbofuran (1563-66-2)	0,010	N
Carbofuran-3-keto (16709-30-1)	0,010	Y
Carboxin (5234-68-4)	0,010	Y
Chinomethionat (2439-01-2)	0,005	N
Chlorbufam (1967-16-4)	0,005	N
Chlordimeform (6164-98-3)	0,005	Y
Chlorfenvinphos (470-90-6)	0,005	Y
Chlorflurenol-Methyl ester (2536-31-4)	0,010	N
Chlorimuron-Ethyl (90982-32-4)	0,010	Y
Chloroxuron (1982-47-4)	0,005	Y
Chlorsulfuron (64902-72-3)	0,010	Y
Chlortoluron (15545-48-9)	0,005	Y
Chromafenozide (143807-66-3)	0,005	Y
Cinosulfuron (94593-91-6)	0,020	Y
Clethodim (99129-21-2)	0,005	N
Clethodim-Sulfone (111031-17-5)	0,005	Y
Clodinafop-Propargyl (105512-06-9)	0,005	Y
Clomazone (81777-89-1)	0,010	Y
Clothianidin (210880-92-5)	0,010	Y
Coumaphos (56-72-4)	0,005	Y
Crimidine (535-89-7)	0,005	Y
Crotoxyphos (7700-17-6)	0,005	Y
Crufomate (299-86-5)	0,005	N
Cyazofamid (120116-88-3)	0,005	Y
Cycloheximide (66-81-9)	0,010	N
Cyclosulfamuron (136849-15-5)	0,010	Y
Cycloxydim (101205-02-1)	0,005	N
Cycluron, OMU (2163-69-1)	0,005	Y
Cyflufenamid: sum of cyflufenamid (Z-isomer) and its E-isomer (180409-60-3)	0,010	Y
Cyflumetofen (400882-07-7)	0,010	Y
Cymoxanil (57966-95-7)	0,010	Y
Cyproconazole A (94361-06-5)	0,005	N
Cyproconazole B (94361-06-5)	0,005	N
Cyprodinil (121552-61-2)	0,005	Y
Cyromazine (66215-27-8)	0,005	N
Daimuron (42609-52-9)	0,010	Y
Demeton-s-methyl (919-86-8)	0,005	N
Demeton-S-Methyl-sulfon (17040-19-6)	0,010	Y
Diazinon (333-41-5)	0,005	Y
Dichlormid (37764-25-3)	0,010	Y
Dichlorprop, 2,4-DP (120-36-5)	0,005	N
Diclobutrazol (75736-33-3)	0,005	Y
Diclofop (40843-25-2)	0,050	N
Diclotophos (141-66-2)	0,005	Y
Diethatyl-Ethyl (38727-55-8)	0,005	N

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food and Feed - Fatty Acids

Detection technique: LC-MS/MS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 12

List nr.: 3

Version nr.: 9

Date issued: 28-05-2022

Valid till: 14-08-2022



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Component/Summation and CASnr	Reporting Limit mg/kg	Accredited (Y/N)
Difenoconazole (119446-68-3)	0,005	Y
Difenofoxuron (14214-32-5)	0,010	Y
Diflubenzuron (35367-38-5)	0,010	Y
Dimefox (115-26-4)	0,005	Y
Dimefuron (34205-21-5)	0,010	Y
Dimethachlor (50563-36-5)	0,010	N
Dimethoate (60-51-5)	0,005	Y
Dimethylphenylformamide, 2,4-, DMF (60397-77-5)	0,010	N
Dimethylphenyl-N-Methylformamidine, N-2,4-, DMPF (33089-74-6)	0,050	N
Dimetilan (644-64-4)	0,005	Y
Dimoxystrobin (149961-52-4)	0,005	Y
Diniconazole (83657-24-3)	0,005	Y
Dioxacarb (6988-21-2)	0,005	Y
Disulfoton-sulfone (2497-06-5)	0,005	Y
Disulfoton-sulfoxide (2497-07-6)	0,010	Y
Dithiopyr (97886-45-8)	0,005	Y
DMS (3984-14-3)	0,005	Y
Dodemorph (1593-77-7)	0,005	Y
EPN (2104-64-5)	0,010	Y
Epoxiconazole (133855-98-8)	0,005	Y
Esprocarb (85785-20-2)	0,005	Y
Etaconazole (60207-93-4)	0,010	Y
Ethaboxam (162650-77-3)	0,005	Y
Ethametsulfuron-Methyl (97780-06-8)	0,005	Y
Ethiofencarb-Sulfoxide (53380-22-6)	0,005	Y
Ethiprole (181587-01-9)	0,005	Y
Ethirimol (23947-60-6)	0,005	Y
Ethoprophos (13194-48-4)	0,005	Y
Ethoxysulfuron (126801-58-9)	0,010	Y
Etrifos (38260-54-7)	0,005	Y
Famphur, Famophos (52-85-7)	0,005	Y
Fenamidone (161326-34-7)	0,005	Y
Fenamiphos (22224-92-6)	0,005	N
Fenamiphos-Sulfoxide (31972-43-7)	0,005	Y
Fenarimol (60168-88-9)	0,010	Y
Fenbuconazole (114369-43-6)	0,005	Y
Fenchlorphos oxon (3983-45-7)	0,005	Y
Fenhexamid (126833-17-8)	0,005	Y
Fenobucarb (3766-81-2)	0,005	Y
Fenoxycarb (72490-01-8)	0,010	Y
Fenpiclonil (74738-17-3)	0,010	Y
Fenpropidin (67306-00-7)	0,005	Y
Fenpropimorph (67564-91-4)	0,005	Y
Fenpyrazamine (473798-59-3)	0,005	Y
Fensulfothion (115-90-2)	0,010	Y
Fensulfothion-Oxon (6552-21-2)	0,005	Y
Fensulfothion-oxon-sulfon (6132-17-8)	0,005	Y
Fensulfothion-sulfone (14255-72-2)	0,010	Y
Fenthion (55-38-9)	0,005	N

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food and Feed - Fatty Acids

Detection technique: LC-MS/MS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 12

List nr.: 3

Version nr.: 9

Date issued: 28-05-2022

Valid till: 14-08-2022



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Component/Summation and CASnr	Reporting Limit mg/kg	Accredited (Y/N)
Fenthion-Oxonsulfone (14086-35-2)	0,005	Y
Fenthion-Oxonsulfoxide (6552-13-2)	0,010	Y
Fenthion-Sulfone (3761-42-0)	0,010	Y
Fenuron (101-42-8)	0,005	Y
Ferimzone (sum of (E),(Z) types) (89269-64-7)	0,005	Y
Fipronil (120068-37-3)	0,005	N
Flonicamid (158062-67-0)	0,005	Y
Flonicamid-Metabolite TFNG (207502-65-6)	0,005	N
Florasulam (145701-23-1)	0,005	Y
Fluacrypyrim (229977-93-9)	0,005	Y
Fluazifop-butyl and Fluazifop-P-butyl (69806-50-4 / 79241-46-6)	0,010	Y
Fluazinam (79622-59-6)	0,005	N
Flubendiamide (272451-65-7)	0,010	N
Flucetosulfuron (412928-75-7)	0,020	Y
Flucycloxuron (113036-88-7)	0,010	Y
Flufenacet (142459-58-3)	0,005	Y
Flumetsulam (98967-40-9)	0,010	Y
Fluometuron (2164-17-2)	0,005	Y
Fluopicolide (239110-15-7)	0,010	N
Fluopyram (658066-35-4)	0,005	Y
Flupyradifurone (951659-40-8)	0,005	Y
Fluquinconazole (136426-54-5)	0,010	Y
Fluridone (59756-60-4)	0,005	Y
Flurochloridone (61213-25-0)	0,010	Y
Flusilazole (85509-19-9)	0,005	Y
Flutriafol (76674-21-0)	0,005	Y
Fluxapyroxad (907204-31-3)	0,005	Y
Fomesafen (72178-02-0)	0,005	N
Forchlorfenuron (68157-60-8)	0,010	Y
Formetanate (22259-30-9)	0,005	Y
Formothion (2540-82-1)	0,005	Y
Fosthiazate (98886-44-3)	0,005	Y
Fuberidazole (3878-19-1)	0,005	Y
Furalaxyl (57646-30-7)	0,005	Y
Furathiocarb (65907-30-4)	0,010	N
Furmecyclox (60568-05-0)	0,005	N
Halofenozide (112226-61-6)	0,010	Y
Haloxypop (69806-34-4)	0,010	Y
Haloxypop-2-ethoxyethyl (87237-48-7)	0,010	Y
Haloxypop-methyl and Haloxypop-p-Methyl (69806-40-2 / 72619-32-0)	0,005	N
Heptenophos (23560-59-0)	0,005	Y
Hexaconazole (79983-71-4)	0,005	Y
Hexaflumuron (86479-06-3)	0,005	N
Hexazinone (51235-04-2)	0,005	Y
Hydramethylnon (67485-29-4)	0,020	Y
Hymexazol, Hydroxyisoxazol (10004-44-1)	0,005	N
Imazalil (35554-44-0)	0,010	Y
Imazamethabenz-Methyl (81405-85-8)	0,010	Y
Imazaquin (81335-37-7)	0,005	Y

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food and Feed - Fatty Acids

Detection technique: LC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 12

List nr.: 3

Version nr.: 9

Date issued: 28-05-2022

Valid till: 14-08-2022



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Component/Summation and CASnr	Reporting Limit mg/kg	Accredited (Y/N)
Imazosulfuron (122548-33-8)	0,010	Y
Imibenconazole (86598-92-7)	0,005	Y
Imidacloprid (138261-41-3)	0,010	Y
Indole-3-Acetic Acid, IAA (87-51-4)	0,005	N
Iprconazole (125225-28-7)	0,010	Y
Iprobenfos (26087-47-8)	0,010	Y
Iprodione (36734-19-7)	0,005	N
Iprovalicarb (140923-17-7)	0,005	Y
Isazofos (42509-80-8)	0,010	Y
Isofenphos (25311-71-1)	0,005	Y
Isofetamid (875915-78-9)	0,010	Y
Isomethiozin (57052-04-7)	0,005	Y
Isoprocarb (2631-40-5)	0,005	Y
Isoprothiolane (50512-35-1)	0,010	Y
Isoproturon (34123-59-6)	0,005	Y
Isopyrazam (881685-58-1)	0,005	Y
Isoxaben (82558-50-7)	0,010	Y
Isoxaflutole (141112-29-0)	0,005	Y
Karbutilate (4849-32-5)	0,005	Y
Kresoxim-methyl (143390-89-0)	0,005	Y
Linuron (330-55-2)	0,005	Y
Malathion (121-75-5)	0,005	N
Mandestrobin (173662-97-0)	0,005	Y
Mandipropamid (374726-62-2)	0,005	Y
MCPA (94-74-6)	0,005	N
MCP, Mechlorprop, Mecoprop (7085-19-0)	0,005	N
Mecarbam (2595-54-2)	0,005	Y
Mefenpyr-diethyl (135590-91-9)	0,005	Y
Mesosulfuron-Methyl (208465-21-8)	0,020	Y
Metalaxyl and Metalaxyl-M (57837-19-1 / 70630-17-0)	0,020	Y
Metalaxyl free acid (87764-37-2)	0,005	N
Metamitron (41394-05-2)	0,005	Y
Metazachlor (67129-08-2)	0,010	Y
Metconazole (125116-23-6)	0,005	Y
Methabenzthiazuron (18691-97-9)	0,005	Y
Methamidophos (10265-92-6)	0,005	Y
Methidathion (950-37-8)	0,005	N
Methiocarb (2032-65-7)	0,005	Y
Methiocarb-sulfon (2179-25-1)	0,010	Y
Methiocarb-sulfoxide (2635-10-1)	0,010	Y
Metholachlor (51218-45-2)	0,005	Y
Methomyl (16752-77-5)	0,010	N
Methoprotrene (841-06-5)	0,005	Y
Methoxyfenozide (161050-58-4)	0,005	Y
Metobromuron (3060-89-7)	0,005	Y
Metominostrobin E (133408-50-1)	0,005	N
Metosulam (139528-85-1)	0,010	Y
Metoxuron (19937-59-8)	0,010	Y
Metrafenone (220899-03-6)	0,010	Y

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food and Feed - Fatty Acids

Detection technique: LC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 12

List nr.: 3

Version nr.: 9

Date issued: 28-05-2022

Valid till: 14-08-2022



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Component/Summation and CASnr	Reporting Limit mg/kg	Accredited (Y/N)
Metribuzin (21087-64-9)	0,005	N
Metsulfuron-Methyl (74223-64-6)	0,010	Y
Mexacarbate (315-18-4)	0,005	Y
Molinate (2212-67-1)	0,010	Y
Monocrotophos (6923-22-4)	0,005	N
Monolinuron (1746-81-2)	0,010	Y
Napropamide (15299-99-7)	0,010	Y
Neburon (555-37-3)	0,005	Y
Nicosulfuron (111991-09-4)	0,010	Y
Nitenpyram (150824-47-8)	0,005	Y
Nitralin (4726-14-1)	0,010	Y
Novaluron (116714-46-6)	0,010	Y
Ofurace (58810-48-3)	0,005	Y
Omethoate (1113-02-6)	0,010	Y
Orbencarb (34622-58-7)	0,005	Y
Orthosulfamuron (213464-77-8)	0,020	Y
Oryzalin (19044-88-3)	0,010	Y
Oxadiazon (19666-30-9)	0,010	Y
Oxadixyl (77732-09-3)	0,020	Y
Oxamyl (23135-22-0)	0,005	N
Oxamyl-Oxime (30558-43-1)	0,010	Y
Oxasulfuron (144651-06-9)	0,010	Y
Oxaziclomefone (153197-14-9)	0,010	Y
Oxycarboxin (5259-88-1)	0,005	Y
Paclobutrazol (76738-62-0)	0,010	Y
Pebulate (1114-71-2)	0,005	Y
Penconazole (66246-88-6)	0,005	Y
Penoxsulam (219714-96-2)	0,010	Y
Pentachlorophenol (87-86-5)	0,005	N
Pethoxamid (106700-29-2)	0,005	Y
Phenthoate (2597-03-7)	0,010	Y
Phorate-Oxon (2600-69-3)	0,005	N
Phorate-Oxonsulfone (2588-06-9)	0,005	Y
Phorate-sulfone (2588-04-7)	0,005	Y
Phorate-sulfoxide (2588-03-6)	0,010	Y
Phosphamidon (13171-21-6)	0,005	Y
Phoxim (14816-18-3)	0,010	Y
Picoxystrobin (117428-22-5)	0,010	Y
Pinoxaden (243973-20-8)	0,005	Y
Pirimicarb (23103-98-2)	0,005	Y
Pirimicarb, Desmethyl- (30614-22-3)	0,005	Y
Pirimiphos-Ethyl (23505-41-1)	0,005	Y
Pirimiphos-Methyl (29232-93-7)	0,005	Y
Pirimiphos-methyl-N-desethyl (67018-59-1)	0,005	Y
Pretilachlor (51218-49-6)	0,050	N
Primisulfuron-Methyl (86209-51-0)	0,010	Y
Probenazole (27605-76-1)	0,005	N
Prochloraz (67747-09-5)	0,005	Y
Prochloraz metabolite BTS44595 (139520-94-8)	0,005	Y

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food and Feed - Fatty Acids

Detection technique: LC-MS/MS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 12

List nr.: 3

Version nr.: 9

Date issued: 28-05-2022

Valid till: 14-08-2022



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Component/Summation and CASnr	Reporting Limit mg/kg	Accredited (Y/N)
Profenofos (41198-08-7)	0,005	Y
Promecarb (2631-37-0)	0,005	Y
Prometryn (7287-19-6)	0,005	Y
Propachlor (1918-16-7)	0,005	Y
Propachlor oxalamic acid (70628-36-3)	0,005	N
Propaquizafop (111479-05-1)	0,010	Y
Propiconazole (60207-90-1)	0,005	Y
Propoxur (114-26-1)	0,005	N
Propyzamide (23950-58-5)	0,005	N
Proquinazid (189278-12-4)	0,005	N
Prosulfuron (94125-34-5)	0,005	Y
Pyracarbolid (24691-76-7)	0,005	Y
Pyraclonil (158353-15-2)	0,005	Y
Pyraclostrobin (175013-18-0)	0,010	Y
Pyrazolynate (58011-68-0)	0,010	Y
Pyrazophos (13457-18-6)	0,010	Y
Pyrazosulfuron-Ethyl (93697-74-6)	0,010	Y
Pyridalyl (179101-81-6)	0,005	N
Pyridate (55512-33-9)	0,010	Y
Pyridate metabolite (40020-01-7)	0,005	N
Pyrifthalid (135186-78-6)	0,010	Y
Pyrimethanil (53112-28-0)	0,005	Y
Pyriminobac-methyl E (147411-69-6)	0,005	Y
Pyriminobac-methyl Z (147411-70-9)	0,005	Y
Pyrimisulfan (221205-90-9)	0,005	N
Pyriofenone (688046-61-9)	0,010	Y
Pyroquilon, 4-Lilolidone (57369-32-1)	0,010	Y
Pyroxasulfone (447399-55-5)	0,005	N
Pyroxsulam (422556-08-9)	0,010	Y
Quinalphos (13593-03-8)	0,010	Y
Quinclorac-methyl ester (84087-33-2)	0,010	N
Rimsulfuron (122931-48-0)	0,020	Y
Rotenone (83-79-4)	0,005	N
Sebuthylazine (7286-69-3)	0,005	Y
Sebuthylazine, Desethyl- (37019-18-4)	0,005	Y
Sethoxydim (74051-80-2)	0,010	Y
Siduron (1982-49-6)	0,005	Y
Simeconazole (149508-90-7)	0,005	Y
Simetryn (1014-70-6)	0,005	Y
Spinetoram A (187166-40-1)	0,005	N
Spinetoram B (187166-40-1)	0,005	N
Spinosyn A (168316-95-8)	0,005	N
Spinosyn D (168316-95-8)	0,005	N
Spirotetramat (203313-25-1)	0,005	N
Spiroxamine (sum) (118134-30-8)	0,005	N
Sulfentrazone (122836-35-5)	0,005	N
Sulfometuron-Methyl (74222-97-2)	0,010	Y
Sulfotep (3689-24-5)	0,005	Y
Sulfoxaflo (946578-00-3)	0,005	N

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food and Feed - Fatty Acids

Detection technique: LC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Annex to declaration of accreditation (scope of accreditation)

Normative document: EN ISO/IEC 17025:2017

Registration number RvA: L 678

Activity nr.: 12

List nr.: 3

Version nr.: 9

Date issued: 28-05-2022

Valid till: 14-08-2022



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Component/Summation and CASnr	Reporting Limit mg/kg	Accredited (Y/N)
Tebufenozide (112410-23-8)	0,010	Y
Tebutam (35256-85-0)	0,010	N
Tebuthiuron (34014-18-1)	0,005	Y
TEPP (107-49-3)	0,005	Y
Terbufos-Sulfone (56070-16-7)	0,005	Y
Terbufos-sulfoxide (10548-10-4)	0,010	N
Terbumeton (33693-04-8)	0,005	Y
Terbutylazin-desethyl (30125-63-4)	0,005	Y
Terbutryn (886-50-0)	0,005	Y
Tetrachlorvinphos (22248-79-9)	0,005	Y
Tetraconazole (112281-77-3)	0,005	Y
Thiabendazole (148-79-8)	0,005	Y
Thiacloprid (111988-49-9)	0,005	Y
Thiazafuron (25366-23-8)	0,010	Y
Thiazopyr (117718-60-2)	0,005	Y
Thidiazuron (51707-55-2)	0,010	Y
Thifensulfuron-methyl (79277-27-3)	0,005	Y
Thifluzamide (130000-40-7)	0,010	Y
Thiobencarb, Benthicarb (28249-77-6)	0,010	Y
Thiocyclam (31895-22-4)	0,005	N
Thiofanox, Thiofanocarb (39196-18-4)	0,005	Y
Thiofanox-Sulfone (39184-59-3)	0,005	Y
Thiofanox-Sulfoxide (39184-27-5)	0,005	Y
Tiadinil (223580-51-6)	0,010	N
Topramezone (210631-68-8)	0,010	N
Tralkoxydim (87820-88-0)	0,005	Y
Triadimefon (43121-43-3)	0,010	Y
Tri-allate (2303-17-5)	0,010	N
Triasulfuron (82097-50-5)	0,010	Y
Triclopyr (55335-06-3)	0,050	N
Tricyclazole (41814-78-2)	0,005	Y
Trietazine (1912-26-1)	0,005	Y
Trifenylin (Fentin) (668-34-8)	0,010	N
Triflumizole (68694-11-1)	0,005	Y
Triflumuron (64628-44-0)	0,005	Y
Triticonazole (131983-72-7)	0,005	Y
Tritosulfuron (142469-14-5)	0,005	Y
Vamidothion (2275-23-2)	0,010	Y
Vernolate (1929-77-7)	0,010	Y