

Current list of Pesticides within the scope of analytical protocol QTI-C-009:**Determination of the content of pesticides and additives****Equivalent to NEN-EN 15662****Scope: Food and Feed****Detection technique: LC-MSMS and/or GC-MSMS****Performance accredited scope in agreement with guideline SANTE/11321/2021**

Normative document: EN ISO/IEC 17025:2017

Laboratory registration number RvA: L 678

QTI List nr.: 2

Version nr.: 15

Date issued: 14-02-2023

Valid till: N/A, current version

**Quality
Testing
Inspection**

The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Component/Summation and CASnr	RL mg/kg	Techn.
^q 1,2,3,4-Tetrachlorobenzene (634-66-2)	0,002	GC
^q 1,2,4,5-Tetrachlorobenzene (95-94-3)	0,002	GC
^q 1-Naphthaleneacetamide (86-86-2)	0,002	GC ¹
^q 2,3,5-Trimethacarb / 3,4,5-Trimethacarb (2655-15-4 / 2686-99-9)	0,002	LC
^q 2,4,5-T (93-76-5)	0,002	LC
^q 2,4,5-T-Methylester (1928-37-6)	0,002	GC
^q 2,4,5-TP, Fenoprop, Fenoprop-butotyl (93-72-1)	0,002	LC
^q 2,4-D (94-75-7)	0,002	LC
^q 2,4-D-2-ethylhexylester (1928-43-4)	0,002	GC
^q 2,4-DB-Methylester (18625-12-2)	0,002	GC ¹
^q 2,4-D-butylester (94-80-4)	0,002	GC ¹
^q 2,4-D-ethylester (533-23-3)	0,002	GC
^q 2,4-D-isobutylester (1713-15-1)	0,002	GC
^q 2,4-D-Methylester (1928-38-7)	0,002	GC
^q 2-Phenylphenol (90-43-7)	0,005	GC
^q 3-Chloroaniline (108-42-9)	0,002	LC
^q 3-Hydroxycarbofuran (16655-82-6)	0,002	LC
^q 4,4-Dibrombenzophenone (3988-03-2)	0,002	GC
^q 4-Chloro-2-methylaniline (95-69-2)	0,002	LC
^q 4-CPA (122-88-3)	0,002	LC
^q 4-Nitrophenol (100-02-7)	0,005	LC
^q Acetamiprid (135410-20-7)	0,005	LC
^q Acetamiprid, N-desmethyl- (190604-92-3)	0,002	LC
^q Acetochlor (34256-82-1)	0,002	GC ¹
Acibenzolar Acid, CGA 210007 (35272-27-6)	0,050	LC
^q Acibenzolar-S-Methyl (135158-54-2)	0,002	LC
^q Aclonifen (74070-46-5)	0,002	GC
Acrinathrin (101007-06-1)	0,002	GC ¹
^q Alachlor (15972-60-8)	0,002	GC ¹
^q Albendazole (54965-21-8)	0,002	LC
^q Aldicarb (116-06-3)	0,002	LC
Aldicarb-Sulfone, Aldoxycarb (1646-88-4)	0,005	LC
^q Aldicarb-Sulfoxide (1646-87-3)	0,002	LC
^q Aldrin (309-00-2)	0,002	GC
^q Ametoctradin (865318-97-4)	0,002	LC
^q Amidosulfuron (120923-37-7)	0,002	LC
Aminocarb (2032-59-9)	0,002	LC
Aminopyralid (150114-71-9)	0,050	LC
^q Amisulbrom (348635-87-0)	0,002	LC
^q Ancymidol (12771-68-5)	0,002	GC
^q Anilofos (64249-01-0)	0,002	LC
^q Anthraquinone (84-65-1)	0,002	GC
^q Aramite (140-57-8)	0,002	GC
^q Aspon (3244-90-4)	0,002	LC
^q Asulam (3337-71-1)	0,002	LC
^q Atraton (1610-17-9)	0,002	LC
^q Atrazine (1912-24-9)	0,002	LCGC
^q Atrazine-desethyl (6190-65-4)	0,002	LC
^q Azamethiphos (35575-96-3)	0,002	GC ¹
^q Azimsulfuron (120162-55-2)	0,002	LC
^q Azinphos-ethyl (2642-71-9)	0,002	LC

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food and Feed

Detection technique: LC-MSMS and/or GC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Normative document: EN ISO/IEC 17025:2017

Laboratory registration number RvA: L 678

QTI List nr.: 2

Version nr.: 15

Date issued: 14-02-2023

Valid till: N/A, current version



**Quality
Testing
Inspection**

The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)
qti-services.com | qti@qti-services.com

Component/Summation and CASnr	RL mg/kg	Techn.
^q Azinphos-methyl (86-50-0)	0,002	LC
^q Aziprotryne (4658-28-0)	0,002	LC
^q Azoxystrobin (131860-33-8)	0,002	GC
^q BAC-10 (965-32-2)	0,002	LC
^q BAC-12 (139-07-1)	0,020	LC
BAC-14 (139-08-2)	0,002	LC
^q BAC-16 (122-18-9)	0,002	LC
^q Benazolin (3813-05-6)	0,002	LC
^q Benazolin-Ethyl (25059-80-7)	0,002	LC
^q Bendiocarb (22781-23-3)	0,005	LC
^q Benfluralin (1861-40-1)	0,002	GC
^q Benfuresate (68505-69-1)	0,002	GC
^q Benodanil (15310-01-7)	0,002	GC
^q Benoxacor (98730-04-2)	0,002	GC ¹
^q Bensulfuron-Methyl (83055-99-6)	0,002	LC
Bentazone (25057-89-0)	0,002	LC
^q Benthiavalicarb-Isopropyl (177406-68-7)	0,002	LC
^q Benthiazole (21564-17-0)	0,002	LC
^q Benzobicyclon (156963-66-5)	0,002	LC
^q Benzovindiflupyr (Solatenol) (1072957-71-1)	0,002	LC
^q Benzoylprop-Ethyl (22212-55-1)	0,002	GC
Bifenazate (149877-41-8)	0,002	LC
^q Bifenthrin (82657-04-3)	0,002	GC
Biphenyl (92-52-4)	0,010	GC
^q Bispyribac (125401-75-4)	0,002	LC
^q Bixafen (581809-46-3)	0,002	LC
^q Boscalid (188425-85-6)	0,002	LCGC
Brodifacoum (56073-10-0)	0,002	LC
^q Bromacil (314-40-9)	0,005	LC
^q Bromfenvinphos-methyl (13104-21-7)	0,002	GC
^q Bromocyclen (1715-40-8)	0,002	GC ¹
^q Bromophos-ethyl (4824-78-6)	0,002	GC
^q Bromophos-methyl (2104-96-3)	0,002	GC
^q Bromopropylate (18181-80-1)	0,002	GC
^q Bromoxynil (1689-84-5)	0,002	LC
^q Bromoxynil-Methyl (3336-39-8)	0,002	GC
^q Bromuconazole a (116255-48-2)	0,002	LC
^q Bromuconazole b (116255-48-2)	0,002	LC
^q Bupirimate (41483-43-6)	0,002	LC
^q Buprofezin (69327-76-0)	0,002	LC
^q Butachlor (23184-66-9)	0,005	LC
^q Butafenacil (134605-64-4)	0,002	LC
^q Butamifos (36335-67-8)	0,002	GC ¹
^q Butocarboxim (34681-10-2)	0,002	LC
Butocarboxim-Sulfoxide (34681-24-8)	0,050	LC
^q Butoxycarboxim (34681-23-7)	0,002	LC
^q Butralin (33629-47-9)	0,002	GC
^q Buturon (3766-60-7)	0,002	LC
^q Butylate (2008-41-5)	0,002	LC
^q Cadusafos (95465-99-9)	0,002	LCGC
^q Carbaryl (63-25-2)	0,002	LC

Current list of Pesticides within the scope of analytical protocol QTI-C-009:**Determination of the content of pesticides and additives****Equivalent to NEN-EN 15662****Scope: Food and Feed****Detection technique: LC-MSMS and/or GC-MSMS****Performance accredited scope in agreement with guideline SANTE/11321/2021**

Normative document: EN ISO/IEC 17025:2017

Laboratory registration number RvA: L 678

QTI List nr.: 2

Version nr.: 15

Date issued: 14-02-2023

Valid till: N/A, current version

**Quality
Testing
Inspection**

The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)
qti-services.com | qti@qti-services.com

Component/Summation and CASnr	RL mg/kg	Techn.
^a Car bendazim (10605-21-7)	0,002	LC
^a Car betamide (16118-49-3)	0,002	LC
Carbofuran (1563-66-2)	0,002	LC
^a Carbofuran-3-keto (16709-30-1)	0,002	LC
^a Carbophenothion (786-19-6)	0,002	LCGC
^a Carbophenothion-Methyl (953-17-3)	0,002	GC
^a Carboxin (5234-68-4)	0,002	LC
^a Carfentrazone-Ethyl (128639-02-1)	0,002	LC
^a Carpropamid (104030-54-8)	0,002	LC
^a Chinomethionat (2439-01-2)	0,002	LC
Chloramben (133-90-4)	0,050	LC
^a Chloramben-methyl ester (7286-84-2)	0,002	GC ¹
^a Chlorantraniliprole (Rynaxapyr) (500008-45-7)	0,002	LC
^a Chlorbenside (103-17-3)	0,002	GC
Chlorbufam (1967-16-4)	0,010	LC
^a Chlordane-cis (5103-71-9)	0,002	GC
^a Chlordane-oxy (27304-13-8)	0,002	GC
^a Chlordane-trans (5103-74-2)	0,002	GC
^a Chlordimeform (6164-98-3)	0,002	LC
^a Chlorethoxyfos (54593-83-8)	0,002	GC
^a Chlorfenapyr (122453-73-0)	0,002	GC
^a Chlorfenprop-Methyl (14437-17-3)	0,002	GC
^a Chlorfenson (80-33-1)	0,002	GC
^a Chlorfenvinphos (470-90-6)	0,002	LC
^a Chlorfluazuron (71422-67-8)	0,002	LC
Chlorflurenol-Methyl ester (2536-31-4)	0,050	LC
Chloridazon (1698-60-8)	0,050	LC
^a Chlorimuron-Ethyl (90982-32-4)	0,002	LC
^a Chlormephos (24934-91-6)	0,002	GC
^a Chlormequat (7003-89-6)	0,002	LC
^a Chlorobenzilate (510-15-6)	0,002	GC ¹
^a Chlorobenzuron (57160-47-1)	0,002	LC
^a Chloroneb (2675-77-6)	0,002	GC
^a Chloropropylate (5836-10-2)	0,002	GC
^a Chloroxuron (1982-47-4)	0,002	LC
^a Chlorpropham (101-21-3)	0,002	GC
^a Chlorpyrifos (2921-88-2)	0,002	LCGC
^a Chlorpyrifos-Methyl (5598-13-0)	0,002	LCGC ¹
^a Chlorsulfuron (64902-72-3)	0,002	LC
^a Chlorthal-Dimethyl (1861-32-1)	0,002	GC
^a Chlorthion (500-28-7)	0,002	GC ¹
^a Chlorthiophos (60238-56-4)	0,002	GC
^a Chlortoluron (15545-48-9)	0,002	LC
^a Chlozolate (84332-86-5)	0,002	GC
^a Chromafenozide (143807-66-3)	0,002	LC
^a Cinidon-Ethyl (142891-20-1)	0,002	LCGC
^a Cinosulfuron (94593-91-6)	0,002	LC
^a Clethodim (99129-21-2)	0,002	LC
^a Clethodim-Sulfone (111031-17-5)	0,002	LC
^a Clodinafop (114420-56-3)	0,005	LC
^a Clodinafop-Propargyl (105512-06-9)	0,002	LC

Current list of Pesticides within the scope of analytical protocol QTI-C-009:**Determination of the content of pesticides and additives****Equivalent to NEN-EN 15662****Scope: Food and Feed****Detection technique: LC-MSMS and/or GC-MSMS****Performance accredited scope in agreement with guideline SANTE/11321/2021**

Normative document: EN ISO/IEC 17025:2017

Laboratory registration number RvA: L 678

QTI List nr.: 2

Version nr.: 15

Date issued: 14-02-2023

Valid till: N/A, current version

**Quality
Testing
Inspection**

The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Component/Summation and CASnr	RL mg/kg	Techn.
Clofentezine (74115-24-5)	0,050	LC
^a Clomazone (81777-89-1)	0,002	LC
Clopyralid (1702-17-6)	0,050	LC
Cloquintocet-Mexyl (99607-70-2)	0,005	GC ¹
Clothianidin (210880-92-5)	0,002	LC ¹
^a Coumaphos (56-72-4)	0,002	LC ¹
^a Crimidine (535-89-7)	0,002	LC
^a Crotoxyphos (7700-17-6)	0,002	LC
^a Crufomate (299-86-5)	0,002	LC
Cyanazine (21725-46-2)	0,002	GC
^a Cyanofenphos (13067-93-1)	0,002	GC
^a Cyanophos (2636-26-2)	0,002	GC
Cyantraniliprole (736994-63-1)	0,002	LC
^a Cyazofamid (120116-88-3)	0,002	LC
^a Cyclanilide (113136-77-9)	0,002	LC
^a Cycloate (1134-23-2)	0,002	GC
^a Cycloheximide (66-81-9)	0,002	LC
^a Cyclosulfamuron (136849-15-5)	0,002	LC
^a Cycloxydim (101205-02-1)	0,002	LC
^a Cycluron, OMU (2163-69-1)	0,005	LC
^a Cyenopyrafen (560121-52-0)	0,002	LC
^a Cyflufenamid: sum of cyflufenamid (Z-isomer) and its E-isomer (180409-60-3)	0,002	LC
^a Cyfluthrin (sum) (68359-37-5)	0,002	GC
^a Cyhalofop-Butyl (122008-85-9)	0,002	GC
^a Cyhalothrin (sum) (Gamma- and Lambda-Cyhalothrin) (76703-62-3/91465-08-6)	0,002	GC
^a Cymoxanil (57966-95-7)	0,002	LC ¹
^a Cypermethrin (sum) (52315-07-8)	0,002	GC
Cyphenothrin (39515-40-7)	0,002	GC
^a Cyproconazole A (94361-06-5)	0,002	LC
^a Cyproconazole B (94361-06-5)	0,002	LC
^a Cyprodinil (121552-61-2)	0,002	LCGC
^a Cyprofuram (69581-33-5)	0,002	GC
Cyromazine (66215-27-8)	0,005	LC
^a Daimuron (42609-52-9)	0,002	LC
^a DDD-OP (53-19-0)	0,002	GC
^a DDD-PP (72-54-8)	0,002	GC
^a DDE-OP (3424-82-6)	0,002	GC
^a DDE-PP (72-55-9)	0,002	GC
^a DDT-OP (789-02-6)	0,002	GC
^a DDT-PP (50-29-3)	0,002	GC
^a DEF/Tribufos (78-48-8)	0,002	GC ¹
^a Deltamethrin (52918-63-5)	0,002	GC
^a Demeton-O (298-03-3)	0,002	LC
^a Demeton-S (126-75-0)	0,002	LC
^a Demeton-s-methyl (919-86-8)	0,005	LC
^a Demeton-S-Methyl-sulfon (17040-19-6)	0,002	LC
^a Desmetryn (1014-69-3)	0,002	GC ¹
^a Di-Allate (2303-16-4)	0,002	LCGC
^a Diazinon (333-41-5)	0,002	LCGC
^a Dicapthon (2463-84-5)	0,002	GC
^a Dichlofenthion (97-17-6)	0,002	GC ¹

Current list of Pesticides within the scope of analytical protocol QTI-C-009:**Determination of the content of pesticides and additives****Equivalent to NEN-EN 15662****Scope: Food and Feed****Detection technique: LC-MSMS and/or GC-MSMS****Performance accredited scope in agreement with guideline SANTE/11321/2021**

Normative document: EN ISO/IEC 17025:2017

Laboratory registration number RvA: L 678

QTI List nr.: 2

Version nr.: 15

Date issued: 14-02-2023

Valid till: N/A, current version

**Quality
Testing
Inspection**

The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)
qti-services.com | qti@qti-services.com

Component/Summation and CASnr	RL mg/kg	Techn.
^a Dichlormid (37764-25-3)	0,002	LC
^a Dichlorobenzamide, 2,6- (2008-58-4)	0,002	GC
^a Dichlorobenzophenone, 2,4- (90-98-2)	0,002	GC ¹
^a Dichlorobenzophenone, 4,4- (90-98-2)	0,002	GC ¹
Dichlorprop, 2,4-DP (120-36-5)	0,002	LC
Dichlorvos (62-73-7)	0,002	LCGC
^a Diclobutrazol (75736-33-3)	0,002	LC
^a Diclocymet (139920-32-4)	0,002	GC
Diclofop (40843-25-2)	0,005	LC
^a Diclofop-Methyl (51338-27-3)	0,002	GC
^a Dicloran (99-30-9)	0,002	GC
^a Dicrotophos (141-66-2)	0,002	LC
^a Dieldrin (60-57-1)	0,002	GC
^a Diethatyl-Ethyl (38727-55-8)	0,002	LC
^a Diethofencarb (87130-20-9)	0,002	GC
^a Difenconazole (119446-68-3)	0,002	LCGC
^a Difenoxuron (14214-32-5)	0,002	LC
^a Difenzoquat methyl sulfate (43222-48-6)	0,002	LC
^a Diflubenzuron (35367-38-5)	0,002	LC
^a Diflufenican (83164-33-4)	0,002	LCGC
^a Dimefox (115-26-4)	0,002	LC
^a Dimefuron (34205-21-5)	0,002	LC
^a Dimethachlor (50563-36-5)	0,002	LCGC
^a Dimethenamid (87674-68-8)	0,002	LC
^a Dimethipin (55290-64-7)	0,002	GC
^a Dimethoate (60-51-5)	0,002	LC
^a Dimethomorph (110488-70-5)	0,002	LC
^a Dimethomorph a (110488-70-5)	0,002	GC
^a Dimethomorph b (110488-70-5)	0,002	GC
Dimethylanilin, 2,4- (95-68-1)	0,002	LC
^a Dimethylphenylformamide, 2,4-, DMF (60397-77-5)	0,002	LC
^a Dimethylvinphos (Z type) (67628-93-7)	0,002	GC
^a Dimetilan (644-64-4)	0,002	LC
^a Dimoxystrobin (149961-52-4)	0,002	LC
^a Diniconazole (83657-24-3)	0,002	LC
^a Dinitramine (29091-05-2)	0,002	GC
Dinobuton (973-21-7)	0,002	GC
^a Dinoseb methyl ether (6099-79-2)	0,002	GC ¹
^a Dinotefuran (165252-70-0)	0,002	LC
^a Dinoterb, DNTBP (1420-07-1)	0,005	LC
^a Dioxabenzofos (Salithion) (3811-49-2)	0,005	LC
^a Dioxacarb (6988-21-2)	0,002	LC
^a Dioxathion (78-34-2)	0,002	GC
^a Diphenamid (957-51-7)	0,002	GC
^a Diphenylamine (122-39-4)	0,002	GC ¹
^a Dipropetryn (4147-51-7)	0,002	GC ¹
Disulfoton (298-04-4)	0,002	GC ¹
^a Disulfoton-sulfone (2497-06-5)	0,002	LC
^a Disulfoton-sulfoxide (2497-07-6)	0,002	LC
^a Ditalimfos (5131-24-8)	0,002	LC
^a Dithiopyr (97886-45-8)	0,002	LC

Current list of Pesticides within the scope of analytical protocol QTI-C-009:**Determination of the content of pesticides and additives****Equivalent to NEN-EN 15662****Scope: Food and Feed****Detection technique: LC-MSMS and/or GC-MSMS****Performance accredited scope in agreement with guideline SANTE/11321/2021**

Normative document: EN ISO/IEC 17025:2017

Laboratory registration number RvA: L 678

QTI List nr.: 2

Version nr.: 15

Date issued: 14-02-2023

Valid till: N/A, current version

**Quality
Testing
Inspection**

The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Component/Summation and CASnr	RL mg/kg	Techn.
^q Diuron (330-54-1)	0,002	LC
DMS (3984-14-3)	0,002	LC
DMST (66840-71-9)	0,002	GC ¹
^q DNOC (534-52-1)	0,002	LC
^q Dodemorph (1593-77-7)	0,002	LC
^q Edifenphos (17109-49-8)	0,002	GC
^q Endosulfan-alpha (959-98-8)	0,002	GC
^q Endosulfan-beta (33213-65-9)	0,002	GC
^q Endosulfansulfate (1031-07-8)	0,002	GC
^q Endrin (72-20-8)	0,002	GC
^q Endrin, Keto- (53494-70-5)	0,002	GC
^q EPN (2104-64-5)	0,002	LC
^q Epoxiconazole (133855-98-8)	0,002	LC
^q EPTC (759-94-4)	0,002	LCGC
^q Esprocarb (85785-20-2)	0,002	LC
^q Etaconazole (60207-93-4)	0,002	LC
^q Ethaboxam (162650-77-3)	0,005	LC
^q Ethalfuralin (55283-68-6)	0,002	GC
^q Ethametsulfuron-Methyl (97780-06-8)	0,002	LC
^q Ethiofencarb (29973-13-5)	0,002	GC ¹
Ethiofencarb-Sulfone (53380-23-7)	0,002	LC
^q Ethion (563-12-2)	0,002	GC
^q Ethiprole (181587-01-9)	0,002	LC
Ethirimol (23947-60-6)	0,002	LC
^q Ethofumesate (26225-79-6)	0,002	GC ¹
^q Ethoprophos (13194-48-4)	0,002	LC
^q Ethoxysulfuron (126801-58-9)	0,005	LC
^q Etofenprox (80844-07-1)	0,002	GC
^q Etoxazole (153233-91-1)	0,002	GC
^q Etridiazole (2593-15-9)	0,002	GC
^q Etrimfos (38260-54-7)	0,002	LC
ETU (96-45-7)	0,010	LC
^q Famoxadone (131807-57-3)	0,002	GC
^q Famphur, Famophos (52-85-7)	0,002	LC
^q Fenamidone (161326-34-7)	0,002	LC
^q Fenamiphos (22224-92-6)	0,002	LC
Fenamiphos sulfone (31972-44-8)	0,002	GC
^q Fenamiphos-Sulfoxide (31972-43-7)	0,002	LC
^q Fenarimol (60168-88-9)	0,002	LC
^q Fenazaquin (120928-09-8)	0,002	LC
^q Fenbuconazole (114369-43-6)	0,002	LC
^q Fenchlorazole-Ethyl (103112-35-2)	0,002	LC
^q Fenchlorphos (299-84-3)	0,002	GC ¹
^q Fenchlorphos oxon (3983-45-7)	0,002	LC
^q Fenfluthrin (75867-00-4)	0,002	GC ¹
^q Fenfuram (24691-80-3)	0,002	GC
^q Fenhexamid (126833-17-8)	0,002	LC
^q Fenitrothion (122-14-5)	0,002	GC ¹
Fenobucarb (3766-81-2)	0,002	LC
^q Fenoxanil (115852-48-7)	0,005	LC
Fenoxaprop and Fenoxaprop-P (95617-09-7 / 113158-40-0)	0,005	LC

Current list of Pesticides within the scope of analytical protocol QTI-C-009:**Determination of the content of pesticides and additives****Equivalent to NEN-EN 15662****Scope: Food and Feed****Detection technique: LC-MSMS and/or GC-MSMS****Performance accredited scope in agreement with guideline SANTE/11321/2021**

Normative document: EN ISO/IEC 17025:2017

Laboratory registration number RvA: L 678

QTI List nr.: 2

Version nr.: 15

Date issued: 14-02-2023

Valid till: N/A, current version

**Quality
Testing
Inspection**

The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)
qti-services.com | qti@qti-services.com

Component/Summation and CASnr	RL mg/kg	Techn.
^a Fenoxaprop-Ethyl and Fenoxaprop-P-Ethyl (66441-23-4 / 71283-80-2)	0,002	LCGC
^a Fenoxycarb (72490-01-8)	0,002	LC
^a Fenpiclonil (74738-17-3)	0,002	LC
^a Fenpropathrin (39515-41-8)	0,002	GC
^a Fenpropidin (67306-00-7)	0,002	LC
^a Fenpropimorph (67564-91-4)	0,005	LC
Fenpyrazamine (473798-59-3)	0,002	LC
^a Fenpyroximate (134098-61-6/111812-58-9)	0,002	LC
^a Fenson (80-38-6)	0,002	GC
^a Fensulfothion (115-90-2)	0,002	LC
^a Fensulfothion-Oxon (6552-21-2)	0,002	LC
^a Fensulfothion-oxon-sulfon (6132-17-8)	0,002	LC
^a Fensulfothion-sulfone (14255-72-2)	0,002	LC
^a Fenthion (55-38-9)	0,002	GC
^a Fenthion-Oxonsulfone (14086-35-2)	0,002	LC
^a Fenthion-Oxonsulfoxide (6552-13-2)	0,002	LC
^a Fenthion-Sulfone (3761-42-0)	0,002	LC
Fenthion-Sulfoxide (3761-41-9)	0,002	GC
^a Fenuron (101-42-8)	0,002	LC
Fenvalerate/Esfenvalerate (RR,SS,RS,SR) (51630-58-1, 66230-04-4)	0,005	GC
^a Ferimzone (sum of (E),(Z) types) (89269-64-7)	0,002	LC
^a Fipronil (120068-37-3)	0,002	LCGC
^a Fipronil-Desulfinyl (205650-65-3)	0,002	GC
^a Fipronil-Sulfide (120067-83-6)	0,002	GC
^a Fipronil-Sulfone (120068-36-2)	0,002	LCGC
^a Flamprop-M-isopropyl (57973-67-8)	0,002	LCGC
^a Flonicamid (158062-67-0)	0,002	LC
^a Flonicamid-Metabolite TFNG (207502-65-6)	0,005	LC
^a Florasulam (145701-23-1)	0,002	LC
^a Fluacrypyrim (229977-93-9)	0,002	LC
^a Fluazifop and Fluazifop-P (69335-91-7 / 83066-88-0)	0,005	LC
^a Fluazifop-butyl and Fluazifop-P-butyl (69806-50-4 / 79241-46-6)	0,002	LCGC
Fluazinam (79622-59-6)	0,002	LC
^a Fluazuron (86811-58-7)	0,002	LC
^a Flubendiamide (272451-65-7)	0,002	LC
^a Flucetosulfuron (412928-75-7)	0,002	LC
^a Fluchloralin (33245-39-5)	0,002	GC
^a Flucycloxuron (113036-88-7)	0,002	LC
^a Flucythrinate (70124-77-5)	0,002	GC
^a Fludioxonil (131341-86-1)	0,002	GC
^a Fluensulfon (318290-98-1)	0,005	LC
^a Flufenacet (142459-58-3)	0,005	LC
Flufenoxuron (101463-69-8)	0,002	LC
^a Flufenpyr-Ethyl (188489-07-8)	0,002	GC
^a Flumetralin (62924-70-3)	0,002	GC
^a Flumetsulam (98967-40-9)	0,002	LC
^a Flumiclorac-pentyl (87546-18-7)	0,002	LC
^a Flumioxazin (103361-09-7)	0,002	LCGC
^a Fluometuron (2164-17-2)	0,002	LC
^a Fluopicolide (239110-15-7)	0,002	LC
^a Fluopyram (658066-35-4)	0,005	LC

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food and Feed

Detection technique: LC-MSMS and/or GC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Normative document: EN ISO/IEC 17025:2017

Laboratory registration number RvA: L 678

QTI List nr.: 2

Version nr.: 15

Date issued: 14-02-2023

Valid till: N/A, current version



**Quality
Testing
Inspection**

The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)
qti-services.com | qti@qti-services.com

Component/Summation and CASnr	RL mg/kg	Techn.
^a Fluorodifen (15457-05-3)	0,002	GC
^a Fluoroglycofene-Ethyl (77501-90-7)	0,002	LC
^a Fluotrimazole (31251-03-3)	0,002	GC ¹
^a Fluoxastrobin (361377-29-9)	0,002	LC
^a Flupyradifurone (951659-40-8)	0,002	LC
^a Flupyrsulfuron-methyl (144740-53-4)	0,005	LC
^a Fluquinconazole (136426-54-5)	0,002	LC
^a Fluralaner (864731-61-3)	0,002	LC
Flurenol-Butyl (2314-09-2)	0,002	GC
^a Fluridone (59756-60-4)	0,002	LC
^a Flurochloridone (61213-25-0)	0,002	LC
Fluroxypyr (69377-81-7)	0,002	LC
Fluroxypyr-Meptyl (81406-37-3)	0,002	LC
^a Flurprimidol (56425-91-3)	0,002	GC ¹
^a Flurtamone (96525-23-4)	0,002	LCGC
^a Flusilazole (85509-19-9)	0,002	LC
^a Flusulfamide (106917-52-6)	0,002	LC
^a Fluthiacet-Methyl (117337-19-6)	0,002	LC
^a Flutolanil (66332-96-5)	0,002	GC
^a Flutriafol (76674-21-0)	0,002	LC
^a Fluvalinate, tau- (102851-06-9)	0,002	GC
^a Fluxapyroxad (907204-31-3)	0,002	LC
^a Fomesafen (72178-02-0)	0,002	LC
^a Fonofos (944-22-9)	0,002	LC
^a Forchlorfenuron (68157-60-8)	0,002	LC
Formetanate (22259-30-9)	0,050	LC
^a Formothion (2540-82-1)	0,002	LC
^a Fosthiazate (98886-44-3)	0,002	LC
^a Fuberidazole (3878-19-1)	0,002	LC
^a Furalaxyl (57646-30-7)	0,002	LC
Furathiocarb (65907-30-4)	0,002	GC
^a Furmecyclox (60568-05-0)	0,002	LC
^a Genite (97-16-5)	0,002	GC
^a Halfenprox (111872-58-3)	0,002	GC
^a Halofenozide (112226-61-6)	0,002	LC
^a Halosulfuron-methyl (100784-20-1)	0,002	LC
^a Haloxyfop (69806-34-4)	0,002	LC
^a Haloxyfop-2-ethoxyethyl (87237-48-7)	0,002	LCGC
^a Haloxyfop-methyl and Haloxyfop-p-Methyl (69806-40-2 / 72619-32-0)	0,002	LCGC
^a HCH-alpha (319-84-6)	0,002	GC
^a HCH-beta (319-85-7)	0,002	GC
^a HCH-delta (319-86-8)	0,002	GC
^a HCH-gamma-(Lindane) (58-89-9)	0,002	GC
^a Heptachlor (76-44-8)	0,002	GC ¹
^a Heptachlorepoxyd cis (1024-57-3)	0,002	GC
^a Heptachlorepoxyd trans (28044-83-9)	0,002	GC
^a Heptenophos (23560-59-0)	0,002	LC
^a Hexachlorobenzene (HCB) (118-74-1)	0,005	GC
^a Hexachlorobutadiene (87-68-3)	0,002	GC
^a Hexaconazole (79983-71-4)	0,002	LC
Hexaflumuron (86479-06-3)	0,010	LC

Current list of Pesticides within the scope of analytical protocol QTI-C-009:**Determination of the content of pesticides and additives****Equivalent to NEN-EN 15662****Scope: Food and Feed****Detection technique: LC-MSMS and/or GC-MSMS****Performance accredited scope in agreement with guideline SANTE/11321/2021**

Normative document: EN ISO/IEC 17025:2017

Laboratory registration number RvA: L 678

QTI List nr.: 2

Version nr.: 15

Date issued: 14-02-2023

Valid till: N/A, current version

**Quality
Testing
Inspection**

The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)
qti-services.com | qti@qti-services.com

Component/Summation and CASnr	RL mg/kg	Techn.
^q Hexazinone (51235-04-2)	0,002	LC
^q Hexythiazox (78587-05-0)	0,002	LC
^q Hydramethylnon (67485-29-4)	0,002	LC
Hymexazol, Hydroxyisoxazol (10004-44-1)	0,050	LC
Imazalil (35554-44-0)	0,010	LC
^q Imazamethabenz-Methyl (81405-85-8)	0,002	LC
^q Imazamox (114311-32-9)	0,002	LC
^q Imazapic (104098-48-8)	0,002	LC
^q Imazapyr (81334-34-1)	0,002	LC
^q Imazaquin (81335-37-7)	0,002	LC
Imazethapyr (81335-77-5)	0,002	LC
^q Imazosulfuron (122548-33-8)	0,002	LC
^q Imibenconazole (86598-92-7)	0,002	LC
^q Imidacloprid (138261-41-3)	0,002	LC
^q Indanofan (133220-30-1)	0,002	GC
^q Indaziflam (950782-86-2)	0,002	GC
^q Indole-3-Butyric Acid, IBA (133-32-4)	0,005	LC
^q Indoxacarb (144171-61-9)	0,002	LCGC ¹
^q Iodofenphos (18181-70-9)	0,002	GC
^q Ioxynil (1689-83-4)	0,002	LC
^q Ioxynil octanoate (3861-47-0)	0,002	GC
^q Ipconazole (125225-28-7)	0,002	LC
^q Iprobenfos (26087-47-8)	0,002	LC
^q Iprodione (36734-19-7)	0,002	LC
Iprovalicarb (140923-17-7)	0,050	LC
Iprovalicarb a (140923-17-7)	0,002	GC
Iprovalicarb b (140923-17-7)	0,002	GC
^q Isazofos (42509-80-8)	0,002	LC
^q Isocarbophos (24353-61-5)	0,002	GC ¹
^q Isodrin (465-73-6)	0,002	GC
^q Isofenphos (25311-71-1)	0,002	LC
^q Isofenphos-Methyl (99675-03-3)	0,002	LCGC
^q Isofenphos-Oxon (31120-85-1)	0,002	GC ¹
^q Isofetamid (875915-78-9)	0,002	LC
^q Isomethiozin (57052-04-7)	0,002	LC
^q Isoprocab (2631-40-5)	0,002	LC
^q Isopropalin (33820-53-0)	0,002	LC
^q Isoprothiolane (50512-35-1)	0,002	LC
^q Isoproturon (34123-59-6)	0,002	LC
^q Isopyrazam (881685-58-1)	0,002	LC
Isoxaben (82558-50-7)	0,002	LC
^q Isoxadifen-Ethyl (163520-33-0)	0,002	GC
^q Isoxaflutole (141112-29-0)	0,002	LC
^q Isoxathion (18854-01-8)	0,002	GC
^q Karbutilate (4849-32-5)	0,002	LC
^q Kresoxim-methyl (143390-89-0)	0,002	LC
^q Lenacil (2164-08-1)	0,002	LC
^q Leptophos (21609-90-5)	0,002	GC
^q Linuron (330-55-2)	0,002	LC
^q Malaoxon (1634-78-2)	0,002	LCGC ¹
^q Malathion (121-75-5)	0,002	GC

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food and Feed

Detection technique: LC-MSMS and/or GC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Normative document: EN ISO/IEC 17025:2017

Laboratory registration number RvA: L 678

QTI List nr.: 2

Version nr.: 15

Date issued: 14-02-2023

Valid till: N/A, current version



**Quality
Testing
Inspection**

The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)
qti-services.com | qti@qti-services.com

Component/Summation and CASnr	RL mg/kg	Techn.
^a Mandestrobin (173662-97-0)	0,002	LC
^a Mandipropamid (374726-62-2)	0,002	LC
^a MCPA (94-74-6)	0,002	LC
MCPA methyl ester (2436-73-9)	0,002	GC
^a MCPA, 2-Ethylhexyl- (29450-45-1)	0,002	GC
^a MCPA, Butyl- (19480-43-4)	0,002	GC
^a MCPA-1-butyl ester (1713-12-8)	0,002	GC
^a MCPA-ethyl ester (2698-38-6)	0,002	GC
^a MCPA-thioethyl (25319-90-8)	0,002	GC
^a MCPP, Mechlorprop, Mecoprop (7085-19-0)	0,002	LC
^a Mecarbam (2595-54-2)	0,002	LC
^a Mefenacet (73250-68-7)	0,002	GC ¹
^a Mefenpyr-diethyl (135590-91-9)	0,002	LC
Mefentrifluconazole (1417782-03-6)	0,005	LC
^a Mepanipyrim (110235-47-7)	0,002	LC
^a Mephosfolan (950-10-7)	0,002	GC ¹
^a Mepiquat (15302-91-7)	0,002	LC
^a Mepronil (55814-41-0)	0,002	GC
^a Mesosulfuron-Methyl (208465-21-8)	0,002	LC
^a Mesotrione (104206-82-8)	0,002	LC
Metaflumizone (139968-49-3)	0,002	LC
^a Metalaxyl and Metalaxyl-M (57837-19-1 / 70630-17-0)	0,002	LCGC ¹
Metamitron (41394-05-2)	0,002	LC
^a Metazachlor (67129-08-2)	0,002	LC
^a Metconazole (125116-23-6)	0,002	LC
^a Methabenzthiazuron (18691-97-9)	0,002	LC
^a Methamidophos (10265-92-6)	0,002	LC
^a Methidathion (950-37-8)	0,002	LC
^a Methiocarb (2032-65-7)	0,002	LC
^a Methiocarb-sulfon (2179-25-1)	0,002	LC
^a Methiocarb-sulfoxide (2635-10-1)	0,002	LC
^a Metolachlor (51218-45-2)	0,002	LC
Methomyl (16752-77-5)	0,002	LC
^a Methoprotetryne (841-06-5)	0,002	LC
^a Methoxyfenozone (161050-58-4)	0,002	LC
^a Metobromuron (3060-89-7)	0,002	LC
^a Metominostrobin E (133408-50-1)	0,002	LC
^a Metosulam (139528-85-1)	0,002	LC
Metoxuron (19937-59-8)	0,002	LC
^a Metrafenone (220899-03-6)	0,002	LC
^a Metribuzin (21087-64-9)	0,002	LC
^a Metsulfuron-Methyl (74223-64-6)	0,002	LC
^a Mevinphos (7786-34-7)	0,002	GC
^a Mexacarbate (315-18-4)	0,002	LC
^a MGK-264 (113-48-4)	0,002	LC
^a Mirex (2385-85-5)	0,002	GC
^a Molinate (2212-67-1)	0,002	LC
Monocrotophos (6923-22-4)	0,050	LC
^a Monolinuron (1746-81-2)	0,002	LC
^a Monuron (150-68-5)	0,002	LC
^a Musk keton (81-14-1)	0,002	GC

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food and Feed

Detection technique: LC-MSMS and/or GC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Normative document: EN ISO/IEC 17025:2017

Laboratory registration number RvA: L 678

QTI List nr.: 2

Version nr.: 15

Date issued: 14-02-2023

Valid till: N/A, current version



**Quality
Testing
Inspection**

The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)
qti-services.com | qti@qti-services.com

Component/Summation and CASnr	RL mg/kg	Techn.
^a Myclobutanil (88671-89-0)	0,002	GC
^a Napropamide (15299-99-7)	0,002	LC
^a Neburon (555-37-3)	0,002	LC
^a Nicosulfuron (111991-09-4)	0,002	LC
^a Nitenpyram (150824-47-8)	0,002	LC ¹
^a Nitralin (4726-14-1)	0,002	GC
^a Nitrapyrin (1929-82-4)	0,002	GC ¹
^a Nitrofen (1836-75-5)	0,002	GC
^a Nitrothal-isopropyl (10552-74-6)	0,002	GC
^a Nonachlor, cis- (5103-73-1)	0,002	GC
^a Nonachlor, trans- (39765-80-5)	0,002	GC
Norflurazon (27314-13-2)	0,002	GC
^a Novaluron (116714-46-6)	0,002	LC
^a Ofurace (58810-48-3)	0,005	LC
^a Omethoate (1113-02-6)	0,002	LC
^a Orbencarb (34622-58-7)	0,002	LC
^a Orthosulfamuron (213464-77-8)	0,002	LC
^a Oryzalin (19044-88-3)	0,005	LC
^a Oxadiargyl (39807-15-3)	0,002	LCGC
^a Oxadiazon (19666-30-9)	0,002	LC
Oxadixyl (77732-09-3)	0,010	LC
^a Oxamyl (23135-22-0)	0,002	LC
^a Oxamyl-Oxime (30558-43-1)	0,002	LC
^a Oxasulfuron (144651-06-9)	0,002	LC
^a Oxathiapiprolin (1003318-67-9)	0,002	LC
^a Oxaziclomefone (153197-14-9)	0,002	LC
^a Oxycarboxin (5259-88-1)	0,002	LC
Oxydemeton-Methyl (301-12-2)	0,002	LC
^a Oxyfluorfen (42874-03-3)	0,002	GC
^a Paclobutrazol (76738-62-0)	0,002	LC
^a Paraoxon (311-45-5)	0,002	GC ¹
^a Paraoxon-methyl (950-35-6)	0,002	GC ¹
^a Parathion (56-38-2)	0,002	GC
^a Parathion-methyl (298-00-0)	0,002	GC ¹
^a PCB 101 (37680-73-2)	0,002	GC
^a PCB 118 (31508-00-6)	0,002	GC
^a PCB 138 (35065-28-2)	0,002	GC
^a PCB 153 (35065-27-1)	0,002	GC
^a PCB 180 (35065-29-3)	0,002	GC
^a PCB 28 (7012-37-5)	0,002	GC ¹
^a PCB 52 (35693-99-3)	0,002	GC
^a Pebulate (1114-71-2)	0,002	LC
^a Penconazole (66246-88-6)	0,002	LC
^a Pencycuron (66063-05-6)	0,002	LC
^a Pendimethalin (40487-42-1)	0,002	GC
^a Penoxsulam (219714-96-2)	0,002	LC
^a Pentachloro-aniline (527-20-8)	0,002	GC
^a Pentachloroanisole (1825-21-4)	0,002	GC
^a Pentachlorobenzene (608-93-5)	0,005	GC
Pentachlorophenol (87-86-5)	0,005	LC
^a Pentachlorothioanisole (1825-19-0)	0,002	GC

Current list of Pesticides within the scope of analytical protocol QTI-C-009:

Determination of the content of pesticides and additives

Equivalent to NEN-EN 15662

Scope: Food and Feed

Detection technique: LC-MSMS and/or GC-MSMS

Performance accredited scope in agreement with guideline SANTE/11321/2021

Normative document: EN ISO/IEC 17025:2017

Laboratory registration number RvA: L 678

QTI List nr.: 2

Version nr.: 15

Date issued: 14-02-2023

Valid till: N/A, current version



**Quality
Testing
Inspection**

The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)
qti-services.com | qti@qti-services.com

Component/Summation and CASnr	RL	
	mg/kg	Techn.
^a Pentanochlor (2307-68-8)	0,002	GC ¹
^a Penthioopyrad (183675-82-3)	0,002	LC
^a Permethrin (52645-53-1)	0,002	GC
^a Perthane (72-56-0)	0,002	GC
^a Pethoxamid (106700-29-2)	0,002	LC
Phenothrin (26002-80-2)	0,002	GC ²
^a Phenthoate (2597-03-7)	0,002	LCGC
^a Phorate (298-02-2)	0,002	GC
^a Phorate-Oxon (2600-69-3)	0,002	LC
^a Phorate-Oxonsulfone (2588-06-9)	0,002	LC
^a Phorate-Oxonsulfoxide (2588-05-8)	0,002	LC
^a Phorate-sulfone (2588-04-7)	0,005	LC
^a Phorate-sulfoxide (2588-03-6)	0,002	LC
^a Phosalone (2310-17-0)	0,002	GC
^a Phosfolan (947-02-4)	0,002	GC ¹
Phosmet (732-11-6)	0,002	LC
^a Phosphamidon (13171-21-6)	0,005	LC
^a Phoxim (14816-18-3)	0,002	LC
Picloram (1918-02-1)	0,050	LC
^a Picolinafen (137641-05-5)	0,002	LCGC
^a Picoxystrobin (117428-22-5)	0,002	LC
^a Pinoxaden (243973-20-8)	0,005	LC
^a Piperonyl-butoxide (51-03-6)	0,002	LCGC ¹
^a Piperophos (24151-93-7)	0,002	GC
^a Pirimicarb (23103-98-2)	0,002	LC
^a Pirimicarb, Desmethyl- (30614-22-3)	0,002	LC
Pirimicarb, Desmethylformamido- (27218-04-8)	0,002	GC ¹
^a Pirimiphos-Ethyl (23505-41-1)	0,002	LC
^a Pirimiphos-Methyl (29232-93-7)	0,002	LCGC
^a Pirimiphos-methyl-N-desethyl (67018-59-1)	0,002	LC
^a Plifenate (21757-82-4)	0,002	GC ¹
^a Prallethrin (23031-36-9)	0,002	GC
^a Primisulfuron-Methyl (86209-51-0)	0,002	LC
^a Probenazole (27605-76-1)	0,002	LC
^a Prochloraz (67747-09-5)	0,002	LC
^a Prochloraz metabolite BTS44595 (139520-94-8)	0,002	LC
^a Prochloraz metabolite BTS44596 (139542-32-8)	0,002	LC
^a Procymidone (32809-16-8)	0,002	GC
Prodiamine (29091-21-2)	0,005	LC
^a Profenofos (41198-08-7)	0,002	LCGC
^a Profluralin (26399-36-0)	0,002	GC ¹
^a Promecarb (2631-37-0)	0,002	LC
^a Prometon (1610-18-0)	0,002	GC
^a Prometryn (7287-19-6)	0,002	LC
Propachlor (1918-16-7)	0,005	LC
Propachlor oxalamic acid (70628-36-3)	0,050	LC ¹
Propamocarb (24579-73-5)	0,002	LC
^a Propanil (709-98-8)	0,002	LC
^a Propaphos (7292-16-2)	0,002	GC
^a Propaquizafop (111479-05-1)	0,002	LC
^a Propargite (2312-35-8)	0,002	GC

Current list of Pesticides within the scope of analytical protocol QTI-C-009:**Determination of the content of pesticides and additives****Equivalent to NEN-EN 15662****Scope: Food and Feed****Detection technique: LC-MSMS and/or GC-MSMS****Performance accredited scope in agreement with guideline SANTE/11321/2021**

Normative document: EN ISO/IEC 17025:2017

Laboratory registration number RvA: L 678

QTI List nr.: 2

Version nr.: 15

Date issued: 14-02-2023

Valid till: N/A, current version

**Quality
Testing
Inspection**

The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)
qti-services.com | qti@qti-services.com

Component/Summation and CASnr	RL mg/kg	Techn.
^q Propazine (139-40-2)	0,002	GC
^q Propetamphos (31218-83-4)	0,002	GC
^q Propham (122-42-9)	0,002	GC
^q Propiconazole (60207-90-1)	0,002	LC
Propisochlor (86763-47-5)	0,002	LC
Propoxur (114-26-1)	0,010	LC ¹
^q Propyzamide (23950-58-5)	0,002	LC
^q Proquinazid (189278-12-4)	0,002	LC
^q Prosulfocarb (52888-80-9)	0,002	GC ¹
^q Prosulfuron (94125-34-5)	0,002	LC
Prothioconazole (178928-70-6)	0,005	LC
^q Prothioconazole-Desthio (120983-64-4)	0,002	GC
^q Prothiofos (34643-46-4)	0,002	GC
^q Prothoate (2275-18-5)	0,002	GC ¹
^q Pyracarbolid (24691-76-7)	0,002	LC
^q Pyraclozil (158353-15-2)	0,002	LC
^q Pyraclostrobin (175013-18-0)	0,002	LC
^q Pyraflufen-Ethyl (129630-19-9)	0,002	GC
Pyrazolynate (58011-68-0)	0,005	LC
Pyrazophos (13457-18-6)	0,002	LC
Pyrazosulfuron-Ethyl (93697-74-6)	0,005	LC
Pyrethrins Cinerin I (25402-06-6)	0,002	LC
Pyrethrins Cinerin II (121-20-0)	0,002	LC
Pyrethrins Jasmolin I (4466-14-2)	0,002	LC
Pyrethrins Jasmolin II (1172-63-0)	0,002	LC
Pyrethrins Pyrethrin I (121-21-1)	0,002	LC
Pyrethrins Pyrethrin II (121-29-9)	0,002	LC
^q Pyributicarb (88678-67-5)	0,002	GC
^q Pyridaben (96489-71-3)	0,002	GC
^q Pyridaphenthion (119-12-0)	0,002	LC
^q Pyridate metabolite (40020-01-7)	0,002	LC
Pyrifeno (88283-41-4)	0,050	LC
^q Pyrifthalid (135186-78-6)	0,002	LC
^q Pyrimethanil (53112-28-0)	0,002	LC
^q Pyrimidifen (105779-78-0)	0,002	GC
^q Pyriminobac-methyl E (147411-69-6)	0,002	LC
^q Pyriminobac-methyl Z (147411-70-9)	0,002	LC
^q Pyriofenone (688046-61-9)	0,002	LC
^q Pyriproxyfen (95737-68-1)	0,002	LCGC
^q Pyroquilon, 4-Lilolidone (57369-32-1)	0,002	LC
^q Pyroxasulfone (447399-55-5)	0,002	LC
^q Pyroxulam (422556-08-9)	0,002	LC
^q Quinalphos (13593-03-8)	0,002	LC
^q Quinlorac-methyl ester (84087-33-2)	0,002	LC
Quinmerac (90717-03-6)	0,050	LC
Quinoclamine, ACNQ (2797-51-5)	0,002	LC
^q Quinoxifen (124495-18-7)	0,002	LCGC
^q Quintozene (82-68-8)	0,002	GC
^q Quizalofop and Quizalofop-P (76578-12-6 / 94051-08-8)	0,002	LC
^q Quizalofop-Ethyl (76578-14-8)	0,002	LC
^q Resmethrin (10453-86-8 / 28434-01-7)	0,002	GC

Current list of Pesticides within the scope of analytical protocol QTI-C-009:**Determination of the content of pesticides and additives****Equivalent to NEN-EN 15662****Scope: Food and Feed****Detection technique: LC-MSMS and/or GC-MSMS****Performance accredited scope in agreement with guideline SANTE/11321/2021**

Normative document: EN ISO/IEC 17025:2017

Laboratory registration number RvA: L 678

QTI List nr.: 2

Version nr.: 15

Date issued: 14-02-2023

Valid till: N/A, current version

**Quality
Testing
Inspection**

The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)

qti-services.com | qti@qti-services.com

Component/Summation and CASnr	RL mg/kg	Techn.
^a Rimsulfuron (122931-48-0)	0,002	LC
Rotenone (83-79-4)	0,002	LC
^a Saflufenacil (372137-35-4)	0,002	LC
^a Sebuthylazine (7286-69-3)	0,002	LC
^a Sebuthylazine, Desethyl- (37019-18-4)	0,002	LC
^a Sedaxane (874967-67-6)	0,002	LC
^a Sethoxydim (74051-80-2)	0,002	LC
^a Siduron (1982-49-6)	0,002	LC
^a Silaneophan (Silafluofen) (105024-66-6)	0,002	GC
^a Silthiofam (175217-20-6)	0,002	GC
^a Simeconazole (149508-90-7)	0,002	LC
^a Simetryn (1014-70-6)	0,002	LC
^a Spinetoram A (187166-40-1)	0,002	LC
^a Spinetoram B (187166-40-1)	0,002	LC
^a Spinosyn A (168316-95-8)	0,002	LC
^a Spinosyn D (168316-95-8)	0,002	LC
^a Spirodiclofen (148477-71-8)	0,005	LC
^a Spiromesifen (283594-90-1)	0,002	LC
^a Spirotetramat (203313-25-1)	0,002	LC
^a Spiroxamine (sum) (118134-30-8)	0,002	LC
^a Sulcotrione (99105-77-8)	0,002	LC
Sulfentrazone (122836-35-5)	0,002	LC
^a Sulfometuron-Methyl (74222-97-2)	0,002	LC
^a Sulfosulfuron (141776-32-1)	0,002	LC
^a Sulfotep (3689-24-5)	0,002	LC
Sulfoxaflor (946578-00-3)	0,005	LC
^a Sulprofos (35400-43-2)	0,002	LCGC
^a Swep (1918-18-9)	0,002	GC
^a Tebuconazole (107534-96-3)	0,002	LC
^a Tebufenozide (112410-23-8)	0,002	LC
^a Tebufenpyrad (119168-77-3)	0,002	LC
^a Tebupirimfos (96182-53-5)	0,002	GC
^a Tebutam (35256-85-0)	0,002	LC
^a Tebutiuron (34014-18-1)	0,002	LC
^a Tecnazene (117-18-0)	0,002	GC
Teflubenzuron (83121-18-0)	0,002	LC
^a Tefluthrin (79538-32-2)	0,002	GC
^a TEPP (107-49-3)	0,002	LC
^a Tepraloxym (149979-41-9)	0,002	LC
Terbacil (5902-51-2)	0,002	GC
^a Terbufos (13071-79-9)	0,002	GC
^a Terbufos-Sulfone (56070-16-7)	0,002	LC
Terbufos-sulfoxide (10548-10-4)	0,002	LC
^a Terbumeton (33693-04-8)	0,002	LC
^a Terbumeton-desethyl (30125-64-5)	0,002	GC ¹
^a Terbutylazin-desethyl (30125-63-4)	0,002	LC
^a Terbutylazine (5915-41-3)	0,002	GC
^a Terbutryn (886-50-0)	0,002	LC
^a Tetrachlorvinphos (22248-79-9)	0,002	LCGC
^a Tetraconazole (112281-77-3)	0,002	LC
^a Tetradifon (116-29-0)	0,002	GC

Current list of Pesticides within the scope of analytical protocol QTI-C-009:**Determination of the content of pesticides and additives****Equivalent to NEN-EN 15662****Scope: Food and Feed****Detection technique: LC-MSMS and/or GC-MSMS****Performance accredited scope in agreement with guideline SANTE/11321/2021**

Normative document: EN ISO/IEC 17025:2017

Laboratory registration number RvA: L 678

QTI List nr.: 2

Version nr.: 15

Date issued: 14-02-2023

Valid till: N/A, current version

**Quality
Testing
Inspection**

The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)
qti-services.com | qti@qti-services.com

Component/Summation and CASnr	RL mg/kg	Techn.
^a Tetramethrin (7696-12-0)	0,002	LC
^a Tetrasul (2227-13-6)	0,002	GC
TFNA (158063-66-2)	0,002	LC
^a Thiabendazole (148-79-8)	0,002	LC
Thiabendazole, 5-Hydroxy- (948-71-0)	0,002	LC
^a Thiacloprid (111988-49-9)	0,002	LC
^a Thiamethoxam (153719-23-4)	0,002	LC
^a Thiazafurion (25366-23-8)	0,002	LC
^a Thiazopir (117718-60-2)	0,002	LC
^a Thidiazuron (51707-55-2)	0,002	LC
^a Thifensulfuron-methyl (79277-27-3)	0,002	LC
^a Thifluzamide (130000-40-7)	0,002	LC
^a Thiobencarb, Benthicarb (28249-77-6)	0,002	LC
Thiocyclam (31895-22-4)	0,050	LC
^a Thiofanox, Thiofanocarb (39196-18-4)	0,002	LC
Thiofanox-Sulfone (39184-59-3)	0,002	LC
^a Thionazin (297-97-2)	0,002	GC
^a Thiophanate-Ethyl (23564-06-9)	0,005	LC
Thiophanate-Methyl (23564-05-8)	0,010	LC
^a Tiadinil (223580-51-6)	0,002	LC
^a Tiocarbazil (36756-79-3)	0,002	GC
^a Tolclofos-Methyl (57018-04-9)	0,002	LCGC ¹
^a Tolfenpyrad (129558-76-5)	0,002	LCGC ¹
Tolyfluanid (731-27-1)	0,050	LC
Topramezone (210631-68-8)	0,050	LC
^a Tralkoxydim (87820-88-0)	0,002	LC
^a Transfluthrin (118712-89-3)	0,002	GC ¹
^a Triadimefon (43121-43-3)	0,002	LCGC
^a Triadimenol (55219-65-3)	0,002	LC
^a Tri-allate (2303-17-5)	0,002	GC
^a Triamiphos (1031-47-6)	0,002	GC
^a Triasulfuron (82097-50-5)	0,002	LC
^a Triazamate (112143-82-5)	0,002	GC
^a Triazophos (24017-47-8)	0,002	GC
^a Tribenuron-Methyl (101200-48-0)	0,002	LC
Trichlorfon (52-68-6)	0,002	LC
^a Trichloronat (327-98-0)	0,002	GC
^a Triclopyr-2-butoxyethyl ester (64700-56-7)	0,002	LC
^a Tricyclazole (41814-78-2)	0,002	LC
^a Tridemorph (24602-86-6)	0,002	LC
^a Tridiphane (58138-08-2)	0,002	GC ¹
^a Trietazine (1912-26-1)	0,002	LC
^a Trifenylin (Fentin) (668-34-8)	0,005	LC
^a Trifloxystrobin (141517-21-7)	0,002	LC
^a Trifloxysulfuron (145099-21-4)	0,002	LC
^a Triflumizole (68694-11-1)	0,002	LC
^a Triflumuron (64628-44-0)	0,002	LC
^a Trifluralin (1582-09-8)	0,002	GC
^a Triflusulfuron-Methyl (126535-15-7)	0,002	LC
Triforin (26644-46-2)	0,050	LC
Trinexapac-ethyl (95266-40-3)	0,050	LC

Current list of Pesticides within the scope of analytical protocol QTI-C-009:
Determination of the content of pesticides and additives
Equivalent to NEN-EN 15662
Scope: Food and Feed
Detection technique: LC-MSMS and/or GC-MSMS
Performance accredited scope in agreement with guideline SANTE/11321/2021
Normative document: EN ISO/IEC 17025:2017
Laboratory registration number RvA: L 678
QTI List nr.: 2
Version nr.: 15
Date issued: 14-02-2023
Valid till: N/A, current version



The next level in food safety

Keenstraat 46 3044CD | Rotterdam (NL)
qti-services.com | qti@qti-services.com

Component/Summation and CASnr	RL	
	mg/kg	Techn.
^Q Triticonazole (131983-72-7)	0,002	LC
^Q Tritosulfuron (142469-14-5)	0,002	LC
^Q Uniconazole (83657-22-1)	0,002	GC ¹
^Q Valifenalate, Valiphenal (283159-90-0)	0,002	LC
^Q Vamidothion (2275-23-2)	0,002	LC
^Q Vernolate (1929-77-7)	0,002	LC
^Q Vinclozolin (50471-44-8)	0,002	GC ¹
^Q Zoxamide (156052-68-5)	0,002	LC

Legend:

Q) within scope of accreditation, RL) Reporting Limit, Numbered references) product(group)s with alternative RLs or excluded product(group)s: 1) Cacao, 2) Rapeseed, 3) Sunflowerseed